
GentleCath Air™ for Women

1/2 with
FeelClean™
Technology

convatec

Convatec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

Manufactured for:
Convatec Inc.
7815 National Service Road, Suite 600
Greensboro, NC 27409 USA
1-800-422-8811

United States of America 1-800-422-8811
France 0800 00 47 67
Brazil +55 11 3529-1800
UK 0800 587 7560
Germany 0800 78 66 200
Italy 800.500.190

© 2024 Convatec Inc.
™ Indicates a Trademark of Unomedical a/s.

Revised: 2024-03
1736691A1

(en) English

GentleCath Air™ for Women

These instructions for use must be read and understood in parallel with the pictorial instructions for use also enclosed prior to using this device. The individual steps in the pictorial instructions, showing the preparation and performance of catheterisation, must be followed.

The device must only be used when prescribed by an appropriately qualified, trained healthcare professional. The patient or carer must follow the advice of the healthcare professional.

Hydrophilic Intermittent Urinary Catheters are intended to provide an intermittent pathway for the drainage of urine from the urinary bladder. The catheter is inserted through the urethra.

Contact a healthcare professional:

- If there is no urine flow or other usage complications occur (e.g., pain, difficult insertion or removal, bleeding, signs of urinary tract infection or discharge)

Indications

- Intermittent drainage of the urinary bladder of females who need assistance with drainage due to conditions causing urinary retention or dysfunction of the urinary system.

User

- Females

Duration of Use

- The catheters have been designed for routine intermittent use. Devices are single use only and must be removed immediately upon drainage of the bladder and disposed of. The overall usage of these devices can be expected to be ongoing repeated multiple times a day.

Environment

- Professional healthcare facilities, homes, work, public bathrooms, out and about

Contraindications

- Do not use in conjunction with irrigation or instillation.

Warnings

- Check the catheter is the correct size as indicated on the packaging.
- Use the catheter immediately upon extraction from the case.
- The catheter is lubricated and ready-to-use once removed from the case. The catheter should never be re-lubricated.
- Do not touch the catheter tube.
- Always slide the catheter slowly and smoothly when inserting and removing. Never force the catheter.
- When removing the catheter, stop if more urine starts to flow. Continue again once urine stops flowing.
- The device is for single use only. Do not re-use or re-insert the device. This may lead to an increased risk of infection, transmission of viruses/bacteria, bleeding, tissue damage or pain.
- Sterility is guaranteed unless the case is damaged or opened prior to use
- Do not use the device if the case is damaged or open prior to use

If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

Storage and Disposal

- Store the product within its case at a recommended room temperature of 10-25 °C (50-77 °F). Keep dry.
- After use, the product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local laws and regulations.

Technical Assistance

If further information or guidance is needed, visit us online at www.gentlecath.com.

Additional paper copies of the Instructions for Use are available from Convatec Limited.

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Cópias em papel adicionais das Instruções de Uso estão disponíveis na Convatec Limited.