# Aquacel™ Ag Surgical

Hydrofiber™ Surgical Cover Dressing with Silver /  
Apósito de Cobertura Quirúrgica de Hydrofiber™  
con Plata / Pansement de recouvrement chirurgical  
Hydrofiber™ avec argent / Penso Cirúrgico de  
Hydrofiber™ com Prata / Hydrofiber™ 含銀手術覆蓋敷料 /  
Hydrofiber™ 含銀手術覆蓋敷料 / Hydrofiber™ Surgical  
커버 드레싱(은 함유)

INSTRUCTIONS FOR USE / INSTRUCCIONES DE USO / MODE  
D'EMPLOI / INSTRUÇÕES DE USO / 使用說明 / 사용 설명서

## en ENGLISH

### Product Description

Aquacel™ Ag Surgical cover dressing is a sterile post operative dressing comprising an inner (wound contact) non woven pad composed of Hydrofiber™ technology and ionic silver sliver-stitch-bonded with nylon and elastane yarns. The pad is held in place between a top layer of Hydrocolloid which is bound to an outer polyurethane film and a windowed skin contact layer consisting of polyurethane film sandwiched between two further layers of Hydrocolloid.

The outer layer of polyurethane film of Aquacel™ Ag Surgical cover dressing provides a bacterial and viral barrier to reduce the risk of infection provided the dressing remains intact and there is no leakage.

The silver in the dressing helps reduce the risk of wound infection. The highly absorbent central Hydrofiber pad absorbs and retains serosanguinous leakage and transforms into a soft gel which maintains a moist environment to support the body's healing process.

The hydrocolloid component of the dressing provides secure, skin friendly adhesion and supports non traumatic removal of the dressing while minimizing damage to the intact surrounding skin.

Aquacel™ Ag Surgical cover dressing is a waterproof, bacterial and viral barrier dressing which is designed to conform to changes in the wound/incision geometry during body movement and post operative edema.

**Indications**  
Under the supervision of a healthcare professional, Aquacel™ Ag Surgical cover dressing may be used for the management of:

Wounds healing by primary intent (e.g. traumatic and elective post operative wounds/incisions) and as an effective barrier to bacterial penetration to help reduce infection.

**Contraindications**  
Aquacel™ Ag Surgical cover dressing should not be used on individuals who are sensitive to the dressing or its components or have had an allergic reaction to the dressing or its components.

**Precautions and Observations**  
**Caution:** Sterility is guaranteed unless the pouch is damaged or opened prior to use.

This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal.

Due to the sterilization process there may be a slight odor on opening the primary packaging.

Do not apply the dressing under tension.

It is not recommended to cut the dressing.

Choose a dressing size to ensure direct contact between the wound/incision and the Hydrofiber pad. The adhesive area should not come into contact with the wound/incision.

Aquacel™ Ag Surgical cover dressing is not designed as a replacement to primary closure methods such as staples or sutures.

Should infection develop during the use of the dressing, appropriate antibiotic therapy should be initiated. The use of Aquacel™ Ag Surgical cover dressing may be continued, but the progress of the wound/incision should be monitored carefully and all treatment should be under medical supervision.

The dressing is not compatible with oil based products such as petrolatum.

While it has been observed that certain silver-containing products may cause skin discoloration following prolonged use, clinical studies with products containing Hydrofiber Technology with ionic silver have shown no such skin discoloration.

**Directions for Use**  
1. Before applying the dressing, cleanse the wound/incision site, rinse well with sterile saline and dry the surrounding skin.

2. Choose a dressing size to ensure direct contact between the wound/incision and the Hydrofiber pad.

3. Remove the dressing from the sterile pack.

4. Initially remove 3/4 of the backing paper along the length of the dressing minimizing finger contact with the central pad and the adhesive surface. Place at one side of the wound/incision line. Gently roll the dressing over the wound/incision line and mold into place.

5. Remove the remaining 1/4 backing paper and gently mold the dressing in place with hand for a secure adhesion. Do not stretch.

6. The translucent hydrocolloid backing allows the clinician to monitor the central Hydrofiber pad and assess when the dressing needs to be changed. The dressing can be left in place for up to 7 days subject to regular clinical assessment and local dressing protocol.

7. All wounds/incisions should be monitored frequently. Remove Aquacel™ Ag Surgical cover dressing when clinically indicated (i.e. leakage, excessive bleeding, suspicion of infection or at 7 days).

**To dress knee incisions. Knee should be approximately at a 30 degree angle for optimal dressing application.**

### Dressing removal

To remove the dressing, press down on the skin with one hand and carefully lift an edge of the dressing with your other hand. Stretch the dressing to break the adhesive seal and remove.

### Storage Precautions

This product is for single use only and is supplied sterile.

If the inner package is damaged, do not use the dressing.

Store at room temperature (10°C - 25°C/50°F - 77°F). Keep dry.

	7-day sustained activity test (simulated wound model)	Zone of inhibition test (spectrum of activity)
Staphylococcus aureus	✓	✓
Pseudomonas aeruginosa	✓	✓
Candida albicans	✓	✓
MRSA		✓
VRE		✓
Serratia marcescens (abtc resistant)		✓
P. aeruginosa (abtc resistant)		✓
Enterobacter cloacae		✓
Klebsiella pneumoniae		✓
Escherichia coli		✓

If further information or guidance is needed, please contact Convatec Professional Services or visit us online at [www.convatec.com](http://www.convatec.com)

Made in Dominican Republic

© 2022 Convatec Inc.

™/® indicates trademarks of Convatec Inc. AQUACEL and HYDROFIBER are registered trademarks of Convatec Inc. in the United States.

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© 2022 ConvaTec Inc.