Kaltostat™ Fortex

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en ENGLISH

PRODUCT DESCRIPTION

Kaltostat™ Fortex Wound Dressing is a soft, white-to-off-white, sterile, non-woven dressing/packing of calcium-sodium alginate fiber. The alginate fibers absorb wound exudate and converts to a firm gel/fiber mat. This gel forms a moist, warm environment at the wound interface. The gel allows trauma-free removal with little or no damage to newly formed tissue. Depending upon the wound, the dressing may be left in place for up to seven days. Kaltostat™ Fortex Wound Dressing is designed to control minor bleeding. This primary dressing should be used with a secondary cover dressing.

INDICATIONS

Kaltostat™ Fortex Wound Dressing is an external wound dressing designed to absorb exudate, control minor bleeding and protect the wound from contamination.

For Over - the - Counter use, Kaltostat™ Fortex Wound Dressing may be used for:

- abrasions
- lacerations
- superficial burns

Under the supervision of a healthcare professional, Kaltostat™ Fortex Wound Dressing may be used for the management of external wounds such as:

- pressure ulcers, venous stasis ulcers, arterial ulcers, diabetic ulcers
- donor sites, post surgical incisions, traumatic wounds.

Kaltostat™ Fortex Wound Dressing may also be used in:

- the local management of bleeding wounds: lacerations, abrasions, nose bleeds, dental extractions, and after surgical wound debridement
- exudate absorption in oncology wounds: fungatingcutaneous tumors, cutaneous metastasis and Kaposi's sarcoma.

CONTRAINDICATIONS

Kaltostat™ Fortex Wound Dressing should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components. Kaltostat™ Fortex Wound Dressing is not indicated for third degree burns or for use as a surgical sponge.

PRECAUTIONS AND OBSERVATIONS

- The use of Kaltostat™ Fortex wound dressing is not recommended for infants less than 12 months of age.
- **Caution:** Sterility is guaranteed unless pouch is damaged or opened prior to use. Single use.
- This wound dressing should not be used with other wound care products without first consulting a healthcare professional.
- Should you observe irritation (reddening, inflammation), maceration (whitening of skin), hypergranulation (excess tissue formation) or sensitivity (allergic reaction), consult a healthcare professional.
- Kaltostat™ Fortex Wound Dressing is designed to create a moist environment that assists in wound healing. If a Kaltostat™ Fortex Wound Dressing has initially formed a gel that is allowed to dry out, removal from the wound can be difficult. Drying out is generally not a problem with exuding wounds. Kaltostat™ Fortex Wound Dressing should be removed with sterile normal saline. Reapplying saline may be necessary to maintain the gel. If the gel dries out, saturate the dried gel with saline to re-hydrate it; this process can take several hours to soften the dried gel.
- The dressing may be used on infected wounds only under the care of a healthcare professional.

In addition, for pressure ulcers, venous stasis ulcers, arterial ulcers, diabetic ulcers, donor sites, post surgical incision, traumatic wounds, oncology wounds and the local management of bleeding wounds:

- Management of the wound types listed above should be under the supervision of a healthcare professional.
- Appropriate supportive measures should be taken where indicated (e.g. use of graduated compression bandaging in the management of venous leg ulcers or pressure relief measures in the management of pressure ulcers).
- Colonization of chronic wounds is common and is not a contraindication to the use of the dressing. The dressing may be used on infected wounds under medical supervision together with appropriate therapy and frequent monitoring of the wound.
- The control of blood glucose, as well as appropriate supportive measures, should be provided with diabetic foot ulcers.
- For oncology wounds, a secondary dressing of high absorbency is recommended.
- Kaltostat™ Fortex Wound Dressing is not intended for use as a surgical sponge. Kaltostat™ Fortex Wound Dressing is not intended to control heavy bleeding. Alternative measures must be considered in those emergency situations where large quantities of blood may be lost.

DIRECTIONS FOR USE

1. Wound Site Preparation and Cleansing:

Before applying Kaltostat™ Fortex Wound Dressing, cleanse the area with an appropriate wound cleanser.

Under the care of a healthcare professional, if necessary, before applying Kaltostat™ Fortex Wound Dressing, the wound should be debrided of excessive necrotic tissue and eschar and irrigated with an appropriate non-toxic cleansing solution.

2. Dressing Preparation and Application

- a. Kaltostat™ Fortex Wound Dressing should be trimmed to the exact size of the wound avoiding any overlap onto the surrounding skin.
- b. For heavily exuding wounds, Kaltostat™ Fortex Wound Dressing should be applied dry onto the wound. An appropriate secondary dressing should be used to secure Kaltostat™ Fortex Wound Dressing in place.
- c. For lightly exuding wounds, Kaltostat™ Fortex Wound Dressing should be placed on the wound and moistened with sterile normal saline. An appropriate secondary dressing should be used to secure Kaltostat™ Fortex Wound Dressing in place.

3. Dressing Change and Removal

- a. On heavily exuding wounds, change the Kaltostat™ Fortex Wound Dressing when strike-through of the secondary dressing occurs or whenever good clinical practice so dictates.
- b. Removal should be trouble free on heavily exuding wounds. Kaltostat™ Fortex Wound Dressing will gel at the wound/dressing interface and lift away.
- c. Removal from lightly exuding wounds may be assisted by saturating the dressing with sterile normal saline.
- d. Cleanse the wound site with an appropriate wound cleanser before applying a new dressing.
- e. As with any dressing, Kaltostat™ Fortex Wound Dressing should be removed from the wound and the wound cleansed at appropriate intervals.

4. Bleeding Wounds

Apply Kaltostat™ Fortex Wound Dressing to the affected area to control minor bleeding. Remove when the bleeding has stopped. **In no instance should the dressing be left in place for more than 7 days.** Apply additional Kaltostat™ Fortex Wound Dressing as outlined above.

Discard any unused portion of the product after dressing the wound.

Protect from light. Store in a cool, dry place.

If further information or guidance is needed, please contact Convatec Professional Services.

Made in UK

© 2025 Convatec

™/® indicates a trademark of the Convatec Group

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Importado por: Convatec PERU S.A.C.