Hydrofiber® Dressing with Strengthening Fibre / Hydrofiber® Verband mit Verstärkenden Fasern / Medicazione in Hydrofiber® con Fibra Rinforzante / Pansement Hydrofiber® avec Fibres Renforcées / Apósito de Hydrofiber® con fibra reforzada / Penso Hydrofiber® com Fibra de Reforço / Hydrofiber® Verband mit Verstärkenden Vezeln / Hydrofiber® Förband med Förstärkande Fiber / Kuituvahvistettu Hydrofiber® -haavasiidos / Hydrofiber® Bandage med förstärkade fibrer / Hydrofiber® Bandasje med forsterkende Fiber / Hydrofiber® Umbind med styrkjandi þráður / Ενίσχυσης Hydrofiber® με ενωγυριστά ίνες / مضادة هيدروفايبر خالطة مسحة بالثياب المقوية

INSTRUCTIONS FOR USE / GEBRAUCHSANLEITUNG / ISTRUZIONI PER L'USO / CONSEILS D'UTILISATION / INSTRUCCIONES DE USO / INSTRUÇÕES DE USO / GEBRUIKSIJNSTRUCTIES / BRUKSANVISNING / KÄYTTÖOHJEET / BRUKSANVISNING / BRUKSANVISNING / LEIDBENINGANG VAN NOTKUN / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / تعليمات استعمال

INTENDED USE
Aquecel® Ribbon Dressings have been designed to be used as a primary dressing for cavity wounds. They are intended to be used under the direction of a healthcare professional as per the indications for use.

INTENDED USES
Aquecel® Ribbon Dressings are intended to be used by healthcare professionals, carers and patients under the direction of a Healthcare Professional.

TARGET PATIENT POPULATION
Aquecel® Ribbon Dressings are intended to be used on patients with one of the wound types listed in the Indications for use.

CLINICAL BENEFITS

Aquecel® Ribbon Dressings absorb wound fluid and bacteria, providing a moist wound healing environment, aiding autolytic debridement and removing dead-space between the wound and dressing interface.

Aquecel® Ribbon Dressings are designed to manage exudate levels which may further damage the wound bed and surrounding skin.

Aquecel® Ribbon Dressings provide a barrier to protect the wound bed from contamination.

INDICATIONS
Aquecel® Ribbon Dressings are indicated for:

- Diabetic foot ulcers
- Pressure ulcers/ulcers
- Surgical wounds
- Traumatic wounds

CONTRAINDICATIONS
Aquecel® Ribbon Dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or any of its components.

PRECAUTIONS AND WARNINGS
• This dressing is not intended for use as a surgical sponge.
• Sterility is guaranteed unless pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or open prior to use and dispose of the device according to local regulations.
• Aquecel® Ribbon Dressings are for single-use only and should not be re-used. Re-use may lead to increased risk of infection, cross-contamination or delayed healing.
• This dressing should not be used with other wound care products without first consulting a healthcare professional.

• Aquecel® Ribbon Dressings can be cut, if required prior to use, across the narrower dimension of the ribbon. Cutting along the longer dimension of the ribbon should not be performed.
• Aquecel® Ribbon Dressings are not compatible with petroleum-based products.
• Due to the sterilisation process there may be a slight odour on opening the primary packaging.
• During the body's normal healing process, non-viable tissue is removed from the wound (autolytic debridement), which could initially make the wound appear larger. If/when formed blood vessels may occasionally produce bloodstained wound fluid following removal of the dressing.
• Aquecel® Ribbon Dressings are designed to be MRI Safe.
• A healthcare professional should be consulted if any of the following are observed during dressing changes, including dressing removal: a) moisture retentive dressing (e.g. film), hyper-granulation (excess tissue formation), signs of infection (increased pain, bleeding, warmth/redness of surrounding tissue, wound exudate), or a change in wound colour and/or odour.
• After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local state and federal laws and regulations.
• If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

ASSOCIATED PRODUCTS & THERAPIES
Aquecel® Ribbon Dressings may be used alone or in combination with other wound care products such as; primary dressings, secondary dressing, securement bandages, compression therapy devices or medicated topical treatments. Use of additional products and therapies should be undertaken with the supervision of a Health Care Professional.

CONTACT AND DURATION OF USE
Aquecel® Ribbon Dressings can be worn for up to 7 days, dressings should be changed either if clinically indicated.

DIRECTIONS FOR USE
• Before applying the dressing, cleanse the wound area with an appropriate wound cleanser.
• Aquecel® Ribbon Dressing with Strengthening Fibre should overlap at least 1 cm onto the skin surrounding the wound.
• When using Aquecel® Ribbon Dressing with Strengthening Fibre in deep cavity wounds, leave at least 2.5 cm outside the wound for easy retrieval.
• Only pack deep wounds up to 80%, as Aquecel® Ribbon Dressing with Strengthening Fibre will expand to fill the wound spaces on contact with wound fluid.
• During the dressing change, gently peel and cover fully with a moisture retentive dressing (e.g., DuoDERM® Extra Thin), or other appropriate secondary dressing. See relevant individual package inserts for complete instructions for use.
• All wounds should be inspected frequently.
• Remove the Aquecel® Ribbon Dressing with Strengthening Fibre when clinically indicated (i.e., leakage, excessive bleeding, or suspicion of infection).
• Aquecel® Ribbon Dressing with Strengthening Fibre can be left in place for up to 7 days, where clinically indicated.

FOR DRY WOUNDS
In addition to the directions for use set forth above:

- Place the Aquecel® Ribbon Dressing with Strengthening Fibre on the wound and wet with sterile water or saline over the wound area only.
- Fully cover the dressing with a moisture retentive dressing such as DuoDERM® Extra Thin to avoid drying out of the dressing and subsequent dressing adherence to the wound.

Discard any unused portion of the dressing
If the immediate product packaging is damaged, do not use.
Store in a cool dry place.
For further information or guidance is needed, please contact Convatec Professional Services.

© 2025 Convatec
™/® indicates a trademark of the Convatec Group

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