Flexi-Seal[™] FECAL COLLECTOR

with hydrocolloid adhesive

Fecal collector with hydrocolloid adhesive / Colector fecal con adhesivo hidrocoloide / Collecteur fécal avec adhésif hydrocolloïde / Fäkalkollektor mit Hydrokolloid-Hautschutz / Collettore fecale con adesivo idrocoloide / Fekal oppsamlingspose med hydrokolloid kleber / Fekalpåse med hydrokolloid vidhäftning / Ulosteenkeräyspussi ihoon kiinnittyväällä hydrokolloidi reunuksella / Fækalopsamlingspose med hydrokolloidklæber / Hidrokolloid yapışkanlı fekal kollektör

ConvTec

ENGLISH

PRODUCT DESCRIPTION

The Flexi-Seal[™] Fecal Collector is comprised of a hydrocolloid adhesive with a "fabric-like" tape border, and a tapered pouch with a drainage port and deodorizing filter.

INDICATION
For use on non-ambulatory patients who are incontinent of loose or semi-formed stool.

CONTRAINDICATION
Use on patients with a known sensitivity to the adhesive, collector, or components.

PATIENT POSITIONING

The patient should be placed in a side lying position, either right or left side, with the upper knee up toward the chest to expose the perianal area.

SKIN PREPARATION

Cleanse perianal skin well with warm water. If soaps or skin cleansers are used, rinse well. Make sure there is no paste residue on the skin from a previous Collector.

Carefully clip any perianal hair growth close to the skin with scissors.

The Hydrocolloid adhesive may be used to fill crevices and skin irregularities. Apply a ring of paste around the edge of the pre-cut opening on the Collector, not on the patient's skin and wait approximately one minute before applying the Collector to skin.

COLLECTOR APPLICATION

If necessary, enlarge pre-cut opening on the barrier. The opening should be just slightly larger than the size of the anus to allow for anal dilatation.

Remove the release paper from the skin barrier only; leave release paper on tape border. Fold the skin barrier in half lengthwise to allow for the best application.

DO NOT TOUCH THE ADHESIVE WITH YOUR FINGERS

Separate the patient's buttocks and place the narrowed portion of the skin barrier against the perineum.

Secure the skin barrier around the anus first and then the area between the anus and the perineal area by pressing firmly and molding the skin barrier with the hand.

Fit the Collector skin barrier to the area between the anus and the coccyx by pressing firmly and molding the skin barrier with the hand.

Remove the release paper from the tape border, smooth into place on the top buttock. Smooth tape into place on the remaining buttock. If the tape border impinges on the vagina or scrotum, carefully trim the tape border with a pair of scissors.

Hold firmly in place for approximately one minute to assure a secure fit.

COLLECTOR POSITION

Smooth Collector Pouch material between patient's legs, avoiding wrinkling of material or pressure areas.

COLLECTOR USAGE

For Liquid Stool

When using bedside drainage: Attach the open cap end of the drainage tap to the end of the drainage tube of a bedside drainage bag.

For ambulation: Disconnect bedside drainage and close cap.

For Semi-Formed Stool

Cut off lower end of Collector immediately above the drainage tap and use the enclosed tail closure to close the Collector.

To attach the tail closure, fold the tail of the Collector over the "knife" edge of the closure *only once*. Hold in place. Press firmly on raised bar until it snaps securely closed. To remove closure, separate "knife" edges from bar.

Open and close the tail closure, as needed, to drain the collector. Do not allow the Collector to become more than one-third full.

Cleaning the Collector

The Collector can be rinsed while in place through the bottom with warm tap water until clean.

REMOVAL OF THE COLLECTOR

Gently peel the skin barrier away from the patient's skin while applying light pressure on the adjacent skin with your free hand.

Fold the skin barrier in half so that it adheres to itself and discard in an appropriate waste receptacle.

The Flexi-Seal[™] Fecal Collector can remain in place up to 48 hours, provided it is secure.

This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use.

If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local laws and regulations.

HOW SUPPLIED

Flexi-Seal[™] Fecal Collector

Box of 10 **Order No. 650078**

For more information about this and other ConvTec products and services, call us toll free at: 1-800-422-8811, Monday through Friday, 8:30 AM to 5 PM, Eastern Time. In Canada, please call 1-800-465-6302. Our Professional Services Department will be happy to answer your questions.

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