Aquacel™ Ag

Surgical

Hydrofiber™ Surgical Cover Dressing with Silver / Chirurgische Hydrofiber™ Wundauflage mit Silber / Medicação Hydrofiber™ com Iodo Argentio per a Ferita Cirúrgica / Pansement Hydrofiber™ Chirurgique de recouvrement à l'argent / Apósito de Cobertura Quirúrgica de Hydrofiber™ com Plata / Pense Quirúrgico de Hydrofiber™ com Plata / Hydrofiber™ Chirurgisch Verband met Zilver / Hydrofiber™ förband med silver för postoperativ bandagering / Hopsella kyläystetty Hydrofiber™ sides, postoperatiivisen sidonnan / Hydrofiber™ bandage med sølv til postoperativ bandagering / Hydrofiber™ bandasje med sølv for postoperativ bandagering / Κυρουργικός επίθεσης εμπαλυσής Hydrofiber™ με άργυρο / جراحة البانسة / Cerrahî Yara Otüsü

INSTRUCTIONS FOR USE / GEBRAUCHSANLEITUNG / ISTRUZIONI PER L'USO / NOTICE D'UTILISATION / INSTRUCCIONES DE USO / INSTRUÇÕES DE USO / GEBRUKSAANWIJZING / BRUKSANVISNING / KÄYTTÖOHJE / BRUGSANVISNING / BRUKSANVISNING / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / تعليمات الاستخدام / KULLANMIN TALIMLARI

ENGLISH

PRODUCT DESCRIPTION

Aquacel™ Ag Surgical cover dressing is a sterile post-operative dressing comprised of an inner (wound contact) non-woven pad composed of Hydrofiber™ technology and a 1.2% iodo silver stitch bonded with nylon and elastane yarns. The pad is held in place between a top layer of Hydrocolloid which is bound to an outer polyurethane film and a windokid skin contact layer consisting of polyurethane film sandwiched between two further layers of Hydrocollod.

The outer polyurethane film of Aquacel™ Ag Surgical cover dressing provides a bacterial and viral barrier to reduce the risk of infection provided the dressing remains intact and there is no leakage.

The 1.2% iodo silver in the dressing helps reduce the risk of wound infection.

This product contains porcine gelatine.

INTENDED USE

Aquacel™ Ag Surgical Cover Dressings have been designed to be used as a primary dressing.

Aquacel™ Ag Surgical Cover Dressings may be used with the direction of a healthcare professional, for the management of acute, surgical wounds which are at risk of infection or show signs of infection.

INTENDED USES

Aquacel™ Ag Surgical Cover Dressings are intended to be used by Healthcare Professionals.

TARGET PATIENT POPULATION

Aquacel™ Ag Surgical Cover Dressings are intended to be worn by patients who have undergone surgery which are at risk of infection.

CLINICAL BENEFITS

Aquacel™ Ag Surgical Cover Dressings provide a physical waterproof, bacterial, and viral barrier to protect the wound.

Aquacel™ Ag Surgical Cover Dressings absorb wound fluid and bacteria, providing a moist wound healing environment, whilst contouring to the wound/sutures/stapes and removing dead-space between the wound and dressing interface.

Aquacel™ Ag Surgical Cover Dressings effectively kill bacteria, yeasts and moulds providing protection against surgical skin infection.

Aquacel™ Ag Surgical Cover Dressings are flexible to allow a full range of movement whilst remaining comfortable to the wearer.

INDICATIONS

Aquacel™ Ag Surgical Cover Dressings are indicated for the management of acute, surgical wounds which are at risk of infection or show signs of infection.

CONTRAINDICATIONS

Aquacel™ Ag Surgical Cover Dressings should not be used on individuals who are sensitive to or who had an allergic reaction to the dressing or any of its components including colophony or its derivatives.

PRECAUTIONS AND WARNINGS

- Sterility is guaranteed unless the pouch is damaged or opened prior to use. Do not use if the integrity of the packaging is damaged or opened prior to use and dispose of the device according to local regulations.
- Aquacel™ Ag Surgical Cover Dressings are for single use only and should not be re-used. Reuse may lead to increased risk of infection and cross contamination and delayed healing.
- This dressing should not be used with other wound care products without first consulting a healthcare professional.
- Aquacel™ Ag Surgical Cover Dressings are not compatible with oil-based products.
- Due to the sterilisation process there may be a slight odour on opening the primary packaging.
- Clinicians/Healthcare professionals should be aware that there are very limited data on prolonged and repeated use of silver containing products, particularly in children and neonates.
- A healthcare professional should be consulted if any of the following are observed during dressing changes, irritation (redness, inflammation), signs of infection (increased pain, heat, hyper-granulation (excess tissue formation), signs of infection (increased pain, bleeding, warmth/redness of surrounding tissue, wound exudate), or a change in wound colour and/or dressing.
- Aquacel™ Ag Surgical Cover Dressings are designed to be MRI-Safe.
- Do not apply the dressing under tension.
- Aquacel™ Ag Surgical Cover Dressings are not designed as a replacement to primary closure methods such as staples or sutures.
- Aquacel™ Ag Surgical Cover Dressings are not designed to be cut.
- It has been observed that silver containing products may cause skin discoloration following prolonged use, however this is often temporary and will reverse when silver contact usage is stopped.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local state and federal laws and regulations.
- If using the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.
- Aquacel™ Ag Surgical Cover Dressings are water resistant however; should not be submerged during bathing or swimming.

ASSOCIATED PRODUCT THERAPIES

Aquacel™ Ag Surgical Cover Dressings are not designed to be used in conjunction with any other products or therapies however, clinical judgement should be used to determine appropriate use where casts or splints are being applied.

CONTACT AND DURATION OF USE

Aquacel™ Ag Surgical Cover Dressings can be worn for up to 7 days, dressings should be changed earlier if clinically indicated. The requirement for Aquacel™ Ag Surgical Cover Dressings should be re-assessed after 14 days and alternative wound management considered where appropriate.

DIRECTIONS FOR USE

To dress knee incisions, knee should be approximately at a 30-degree angle for optimal dressing application.

1. Before applying the dressing, cleanse the wound/incision site with an appropriate cleanser, rinse and pat dry with sterile saline and dry the surrounding skin.
2. Create a dressing size to ensure direct contact between the wound/incision and the Hydrofiber™ technology pad.
3. Remove the dressing from the sterile pack
4. Initially remove 3/4 of the backing paper from along the length of the dressing minimising finger contact with the central pad and the adhesive surface. Place at one side of the wound/incision line and gently roll the dressing over the wound/incision line and mould to place.
5. Remove the remaining 1/4 backing paper and gently mould the dressing over the wound/incision line. Mould the dressing in place for a secure adhesion. Do not stretch.
6. The translucent polyurethane and Hydrocollod backed materials allow the clinician to monitor the central Hydrofiber™ technology pad and assess when the dressing needs to be changed. The dressing can be left in place for up to 7 days subject to regular clinical assessment and local dressing protocol.
7. All wound/incisions should be monitored frequently. Remove Aquacel™ Ag Surgical cover dressing when clinically indicated (i.e. leakage, excessive bleeding, suspicion of infection or at 7 days).

Dressing Removal

To remove the dressing, press down on the skin with one hand and carefully lift an edge of the dressing with your other hand. Stretch the dressing to break the adhesive seal and remove.

Store at room temperature (10°C - 25°C/50°F - 77°F). Keep dry.

<https://e.europa.eu/tools/edumad>. This is the Summary of Safety and Clinical Performance (SSCP) location after the launch of European Database on Medical Devices/ EDUMAD.

If further information or guidance is needed, or to report an adverse event, please contact Convatec®. Contact details can be found at www.convatec.com or use the contact details at the end of this document.

Made in Dominican Republic

© 2025 Convatec Inc.
™Indicates a trademark of Convatec Inc.

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™ indica un marchio di fabbrica di Convatec Inc.