Detachable ostomy system consisting of a pouch and a stoma appliance. The system is made of high quality materials and is designed for long-term use.

Two Piece Modular Ostomy System

INSTRUCTIONS FOR USE
(See also use with the same size NATURA® pouch)

INTENDED USE
Management of stoma output.

INDICATIONS FOR USE
For use with a colostomy, ileostomy or urostomy.

PRECAUTIONS AND OBSERVATIONS
Caution: Patients with a flange ostomy, when associated with adhesives and the absence of a physical reservoir, should know that there are reported cases of adhesive failure in the urinary tract. Initial appropriate clinical management and education are advised. The continued use of the adhesive product if a mass is detected.

Important: A medical professional should assess the need for skin barriers with convalescent prior to use.

Directions for use:

1. Skin Preparation
Your hands and the skin surrounding the stoma must be kept free from any oils substances before applying.

2. Skin Barrier Application
Do not remove the clear covering or release liner from the skin barrier. Skin barrier should be released to ensure proper adhesion to the stoma.

3. Personalize the skin barrier opening, place your thumb(s) finger(s) at the edge of the stoma opening to gently roll the adhesive back. You may thumb(s) finger(s) to continue to roll the adhesive back until you create an opening of the desired size and shape of your stoma. The adhesive must be rolled, do not stretch, pull or cut the adhesive to make it fit. Apply immediately following molding.

Notes: 1) The tape collar may cause skin irritation for users wearing digastric appliances. If irritation occurs, cut the collar away from the skin barrier.

2) Protect the stoma from trauma or injury.

3) Applying the Skin Barrier from the Side
Remove release liner from skin barrier. Center the opening of the skin barrier over the stoma. Apply skin barrier against the skin around the stoma for 30 seconds. If a flange border is present, remove remaining release liner. Smooth in place and remove wrinkles. Press the rolled adhesive to the stoma down in a circular motion for 30 seconds.

Wear as long as comfortable and secure. Wear time may vary.

4) Applying the Skin Barrier from the Front
Separate pouch flims to allow it to breathe. Line up the skin barrier, top and pouch flange and press on the bottom of the flange. Then move the top of the flange, pressing as you go around. Flip top pouch to ensure secure connection.

5) Closing the Pouch
A) Detachable Pouches with INVISOLENE® Tail Closure (Figure 4)
Fold tail closure beginning with the widest section and ending in the narrowest. Press interlocking closures together across the width of the tail and separate until you feel them close. If desired, turn pouch inside-out to conceal tail closure.

B) Detachable Pouches with Tail Cap (Figure 9)
C) Detachable Pouches with Fold Over Top (Figure 11)

6) Emptying the Pouch
A) Detachable Pouches with INVISOLENE® Tail Closure (Figure 4)
B) Detachable Pouches with Tail Cap (Figure 9)
C) Detachable Pouches with Fold Over Top (Figure 11)

7) Connection to Inlet Drainage System
ACUSAFE™ Top (Figure 12)
Fold Over Top (Figure 13)

8) Information About Pouches with Filters
Close the filter cover before cleaning, bathing, swimming, etc. Remove after water activity.

9) Removing The Pouch (Figure 14)
10) Removing The Skin Barrier (Figure 15)

11) Disposal
This device is for single-use only and should not be reused. Do not reuse any leak or infection risk or other contamination. Physiological fluids in the device may no longer be optimal for reuse.

If using the use of this device as a result of its use a serious incident has occurred, please report to the manufacturer and your national authority.

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local laws and regulations.

© 2025 ConvaTec Inc. All rights reserved.
** indicates a trademark of ConvaTec Inc.

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