# AQUACEL™ Foam

**INSTRUCTIONS FOR USE / GEBRAUCHSANWEISUNG / ISTRUZIONI PER L'USO / NOTICE D'UTILISATION / INSTRUCCIONES DE USO / INSTRUÇÕES PARA UTILIZAÇÃO / GEBRUIKSINSTRUCTIES / BRUKSANVISNING / KÄYTTÖOHJEET / BRUKSANVISNING / BRUKSANVISNING / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / تعليمات الاستخدام**

**Adhesive and Non Adhesive Hydrofoam™ Foam Dressing / Adhäsiv/nicht-adhäsiv Hydrofober™ Schaumverband / Medicazione Hydrofober™ adesiva e non adesiva in schiuma / Pansement hydrocellulaire, issu de la technologie Hydrofober™ adhésif et non adhésif / Espuma de Adhêsivo e Não Adhêsivo / Adhesivo e no Adhêsivo / Preso de espuma de Tecnologia Hydrofober™ aderente e não aderente / Adhesiv en nicht-adhësiv Hydrofoam™ schuimverband / Vidhäftande och icke vidhäftande skumförband med Hydrofober™ teknologi / Kiinnittävä ja kiinnittämättömät Hydrofober™ vaahdotusohjeet / Klebende og ikke klebende skumbandage med Hydrofober™ teknologi / Klebende og ikke-klebende Hydrofober™ skumbandage / Αυτοκόλλητο και μη-αυτοκόλλητο Αποκλειστικό Επικάλυψη με Hydrofober™ και μη-αυτοκόλλητο Αποκλειστικό Επικάλυψη με Hydrofober™**

## ENGLISH

**PRODUCT DESCRIPTION**
AQUACEL™ Foam dressing is a sterile, Hydrofober™ foam wound dressing, consisting of:
• a waterproof, breathable outer polyurethane film which provides bacterial and viral barrier
• a soft, multi-layered, absorbent central pad containing polyurethane foam and Hydrofober™ (Sodium carboxymethylcellulose) wound contact layer. The pad absorbs and retains moisture, exudates and bacteria
• AQUACEL™ Foam Adhesive has a silicone adhesive skin contact border to adhere the dressing to the patient.

**CLINICAL BENEFITS**
AQUACEL™ Foam is designed to protect fragile skin from moisture, shear and friction damage. AQUACEL™ Foam is designed to provide a moist wound healing environment, aiding autolytic debridement and removing dead space between the wound and dressing interface. AQUACEL™ Foam is designed to manage excess exudate levels which may further damage the wound bed and surrounding skin.

**INTENDED USE**
AQUACEL™ Foam dressings may be used, with the consultation of a healthcare professional, for the management of existing non-healing wounds and protection of intact skin. AQUACEL™ Foam dressing may be included in a comprehensive protocol of care to protect intact skin against breakdown. AQUACEL™ Foam dressings may be used as a primary or secondary dressing.

**INDICATIONS**
AQUACEL™ Foam dressing is indicated for the treatment of:
• Leg ulcers
• Pressure Ulcers
• Diabetic Foot Ulcers
• Partial Thickness Wounds
• Partial thickness burns
• Traumatic wounds

**INTENDED USERS**
AQUACEL™ Foam dressings are to be used by health care professionals, patients and carers.
**TARGET PATIENT POPULATION**
AQUACEL™ Foam dressings are intended to be used on patients with one of the wound types listed in the indications or those who require skin protection from breakdown.

**CONTRAINDICATIONS**
AQUACEL™ Foam dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or components stated above.

**PRECAUTIONS AND OBSERVATIONS**
Sterility is guaranteed unless the pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or opened prior to use and dispose of the device according to local regulations. Do not use after stated expiry date.

AQUACEL™ Foam dressings are for single use only, do not reuse. Reuse may lead to increased risk of infection and cross-contamination and delayed healing. For highly exuding wounds, AQUACEL™ Foam may be used as a secondary dressing in conjunction with a primary dressing, such as AQUACEL™ Extra. This dressing should not be used with other wound care products without first consulting a healthcare professional.
Adhesive dressing is cut, additional tape may be required to secure in place.
Additional tape, or other means of fixation, will be required to secure the AQUACEL™ Foam Non-Adhesive dressing in place.
AQUACEL™ Foam dressings are water resistant however, should not be submerged during normal swimming.
AQUACEL™ Foam dressings are not compatible with oil-based products.
AQUACEL™ Foam dressings are not compatible with oxidising agents such as hydrogen peroxide or hypochlorite solutions.
During the body's normal healing process, non-viable tissue is removed from the wound (debridement), which could initially make the wound appear larger.
A healthcare professional should be consulted if any of the following are observed during dressing changes, irritation (reddening, inflammation), maceration (whitening of skin), hypergranulation (excess tissue formation), signs of infection (increased pain, bleeding, warmth) or signs of surrounding tissue damage, which could indicate that the wound is not healing.
After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local state and federal laws and regulations.
If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**CONTACT AND DURATION OF USE**
AQUACEL™ Foam dressings can be worn for up to 7 days, dressings should be changed earlier if clinically indicated, and are safe for continuous use.

**DIRECTIONS FOR USE**
**1. Storage**
Store at room temperature (10° C - 25° C/50° F - 77° F). Protect from light. Keep dry.

**2. Wound Site Preparation and Skin Cleaning:**
Before applying the dressing, cleanse the area with an appropriate wound cleanser and dry the surrounding skin.

**3. Dressing Preparation and Application:**
a. When applying to a wound choose a dressing size to ensure that the central absorbent pad covers the wound (1½ inch) beyond the wound.
b. When applying to intact skin choose a dressing which covers the entire vulnerable area.
c. Remove the dressing from the sterile pack.
d. For AQUACEL™ Foam Adhesive only-Remove the release liner, minimising finger contact with the adhesive surface.
e. Hold the dressing over the wound or area of skin to be protected and line up the centre of the dressing with the centre of the wound or area of skin to be protected.
f. Ensure the dressing is not stretched upon application. For AQUACEL™ Foam Adhesive smooth down the adhesive border. Unexposed joints are sufficiently fixed to allow patient mobility.
g. The dressing may be cut, additional tape or film dressing should be used to create a seal and ensure a bacterial barrier.

**4. Dressing Removal:**
a. The dressing should be changed when clinically indicated (the dressing should not be allowed to get over saturated) or if it becomes damaged/soiled.
b. If removing the dressing is difficult, the dressing should be saturated with sterile saline or water and removed slowly.
c. To remove the dressing, press down gently on the skin and carefully lift one corner of the dressing, continue until all edges are free. Carefully lift away the dressing and discard as per local clinical practices.
d. Only during preventative care, dressing corner can be lifted for inspection of skin and resealed. Dressing must be changed if skin is broken.

Visit us online at [www.convatec.com](http://www.convatec.com).
Made in UK
© 2025 Convatec Inc.
\*\*TM indicates a trademark of Convatec Inc. AQUACEL, the AQUACEL logo and Hydrofober are registered trademarks

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