# AQUACEL<sup>®</sup> Pro Foam Pansement Mousse

## INSTRUCTIONS FOR USE / INSTRUCCIONES DE USO / MODE D'EMPLOI

### AQUACEL<sup>®</sup> Foam Pro Adhesive Silicone Hydrofiber<sup>™</sup> Foam dressing / AQUACEL<sup>®</sup> Foam Pro Apósito de espuma de Hydrofiber<sup>™</sup> con adhesive de silcona / Pansement mousse Hydrofiber<sup>MC</sup> avec adhésif en silicone AQUACEL<sup>MD</sup> Foam Pro

en ENGLISH

#### PRODUCT DESCRIPTION

AQUACEL<sup>®</sup> Foam Pro dressing is a sterile Hydrofiber<sup>™</sup> foam wound dressing consisting of a waterproof outer polyurethane film, a multi-layered absorbent pad, and silicone adhesive.

The multi-layered absorbent dressings conform to the skin and wound and contain a layer of polyurethane foam, a non-woven layer of Hydrofiber<sup>™</sup> (sodium carboxymethylcellulose) and a perforated silicone adhesive wound contact layer.

The outer film layer provides a waterproof viral and bacterial barrier which helps protect the wound from external contaminants, reducing the risk of infection.

The film also helps to manage the moisture vapor transmission of the exudate absorbed by the dressing.

The foam and Hydrofiber<sup>™</sup> materials within the dressing absorb high amounts of wound fluid and bacteria. The Hydrofiber<sup>™</sup> layer helps maintain a moist environment in the wound (which supports the body's healing process) and aids in the removal of non-viable tissue from the wound (autolytic debridement) without damaging viable tissue.

The dressing's silicone adhesive wound contact layer provides secure, skin friendly gentle adhesion.

AQUACEL<sup>®</sup> Foam Pro dressing may be used as a secondary dressing. It may be used alone or in combination with other wound care products and can be cut to shape and size, as directed by a healthcare professional. If the dressing is cut, additional tape may be required to secure it in place.

AQUACEL<sup>®</sup> Foam Pro adhesive dressing film backing acts as a barrier to the wound against bacterial and blood-borne viral pathogens (e.g. HIV and Hepatitis viruses).

The use of this dressing neither guarantees nor warrants against AIDS or Hepatitis virus transmission.

#### INDICATIONS

AQUACEL<sup>®</sup> Foam Pro dressings may be used for the management of both chronic and acute wounds, such as:

- Leg ulcers, pressure ulcers (Stage II-IV) and diabetic ulcers;
- Surgical wounds (post-operative, donor sites, dermatological);
- Partial thickness (second degree) burns;
- Traumatic or surgical wounds left to heal by secondary intention such as dehisced surgical incisions;
- Surgical wounds that heal by primary intent, such as dermatological and surgical incisions (e.g. orthopedic and vascular);
- Local management of wounds that are prone to bleeding, such as wounds that have been mechanically or surgically debrided and donor sites;
- Abrasions;
- Lacerations;
- Minor cuts;
- AQUACEL<sup>®</sup> Foam Pro dressing may be included in a comprehensive protocol of care to protect against skin breakdown.

#### CONTRAINDICATIONS

AQUACEL<sup>®</sup> Foam Pro dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.

#### PRECAUTIONS AND OBSERVATIONS

**Caution:** Sterility is guaranteed unless pouch is damaged or opened prior to use.

This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use. This wound dressing should not be used with other wound care products without first consulting a healthcare professional.

During the body's normal healing process, non-viable tissue may be removed (autolytic debridement) from the wound, which could make the wound appear larger after the first few dressing changes.

Should you observe irritation (reddening, inflammation), maceration (whitening of skin), or hypergranulation (excess tissue formation), consult a healthcare professional.

The wound should be inspected during dressing changes for (1) signs of infection (increased pain, bleeding, warmth/redness of surrounding tissue, wound exudate), (2) there is a change in wound color and/or odor, (3) any other unexpected symptoms occur (e.g. maceration or hypergranulation).

**The use of AQUACEL<sup>®</sup> Foam Pro dressing has not been studied in wounds due to herpes simplex or impetigo.**

**In addition, for leg ulcers, pressure ulcers, diabetic ulcers, second degree burns, donor sites, surgical or traumatic wounds left to heal by primary or secondary intention:**

- Treatment of the wound types listed above should only be under the supervision of a healthcare professional.
- Appropriate supportive measures should be taken where indicated (e.g. use of graduated compression bandaging in the management of venous leg ulcers or pressure relief measures in the management of pressure ulcers).
- Colonization of chronic wounds is common and is not a contraindication to the use of the dressing. The dressing may be used on infected wounds under medical supervision together with appropriate therapy and frequent monitoring of the wound.
- The control of blood glucose, as well as appropriate supportive measures, should be provided with diabetic foot ulcers.

AQUACEL<sup>®</sup> Foam Pro dressings are not compatible with oil based products or emollients such as petrolatum.

Do not use AQUACEL<sup>®</sup> Foam Pro dressings with oxidising agents such as hydrogen peroxide or hypochlorite solutions.

The use of dressings as part of a prophylactic therapy does not preclude the need to continue to develop and follow a comprehensive protocol of care to protect against skin breakdown i.e. support surfaces, positioning, nutrition, hydration, skin care and mobility

#### DIRECTIONS FOR USE

**1. If the immediate sterile product pouch is damaged, do not use.**

#### 2. Wound Site Preparation and Skin Cleansing:

Before applying the dressing, cleanse the wound area and skin with an appropriate wound cleanser and dry the surrounding skin. Ensure surrounding skin is free from any emollients or oil based products that might affect dressing adhesion.

#### 3. Dressing Preparation and Application:

- Choose a dressing size to ensure that the central absorbent pad is 1cm (1/3 inch) larger than the wound area.
- Remove the dressing from the sterile pack, minimizing finger contact with the adhesive surface where applicable. Remove the release liner.
- The dressing can be cut to shape for convenience.
- Hold the dressing over the wound or area of skin to be protected and line up the center of the dressing with the center of the wound or area of skin to be protected. Place the pad directly over the wound and smooth down the adhesive border.
- An appropriate retention bandage or tape should be used to secure the dressing in place if required or if the dressing has been cut.

f. Discard any unused portion of the product after dressing the wound.

#### 4. Dressing Removal:

- The dressing should be changed when clinically indicated (i.e.: leakage, bleeding, increased pain, suspicion of infection). Maximum recommended wear time is up to seven days.
- The wound should be cleansed at appropriate intervals.
- To remove the dressing, press down gently on the skin and carefully lift one corner of the dressing. Continue until all edges are free. Carefully lift away the dressing and discard according to local clinical protocols.

Store at room temperature (10° C - 25° C/50° F - 77° F). Protect from light. Keep dry.

If further information or guidance is needed, please contact ConvaTec Professional Services or visit us online at [www.convatec.com](http://www.convatec.com).

Made in UK

© 2023 ConvaTec Inc.

<sup>®/™</sup> indicates a trademark of ConvaTec Inc. AQUACEL, the AQUACEL logo and Hydrofiber are registered trademarks in the United States.

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**ConvaTec Limited**
Importado por: CONVATEC PERU S.A.C.