DuoDERM™ CGF™

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© 2021 ConvaTec Inc.

® indicates a trademark of ConvaTec Inc.

Wound Dressing


For Single Use

Gamma Sterilized

Consult Instructions for Use

Store at room temperature
Keep dry. Avoid refrigeration and exposure to high humidity

Do not use if packaging is damaged.

Lot

Expiry Date

Abbreviation for authorised representative in Switzerland

Importer

PRODUCT DESCRIPTION

DuoDERM CGF is a hydrocolloid dressing consisting of an inner (wound contact) layer of hydrocolloids contained within an adhesive polymer matrix and an outer layer of polyurethane foam. The dressing provides a moist wound environment that favors the healing process by aiding autolytic debridement, facilitating the migration of epithelial cells, and allowing non-traumatic removal of the dressing without damaging newly formed tissue. The dressing may be used as a primary dressing or a secondary dressing in combination with other wound dressings. The adhesive layer contains elastomeric polymers which enhance the dressing's ability to contain wound exudate by forming a cohesive gel.

INDICATIONS

Chronic Wounds - Pressure Ulcers (Stage I - IV), Leg Ulcers. **Acute Wounds** - Traumatic Wounds (Minor Abrasions, Lacerations), Burns (First and Second degree), Dermatological Excisions, Surgical Wounds (Post-operative Wounds, Donor sites).

CONTRAINDICATION


Not intended for use on patients with known sensitivity to the dressing or its components.

PRECAUTIONS AND OBSERVATIONS

This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use.

1) Appropriate supportive measures should be taken where indicated (e.g., use of graduated compression bandaging in the management of venous leg ulcers, or pressure relief measures in the management of pressure ulcers). 2) The autolytic debriding environment produced by the dressing may result in the wound appearing larger after the first few dressing changes. 3) Colonization of chronic wounds is common and is not a contraindication to the use of the dressing. The dressing may be used on infected wounds under medical supervision together with appropriate therapy and frequent monitoring of the wound. In the presence of a clinical anaerobic infection, occlusive therapy is not recommended. 4) Rarely, irritation and/or maceration of the skin, hypergranulation and allergic reactions have been reported. 5) Frequent dressing changes on wounds with compromised peri-ulcer skin is not recommended.

APPLYING DUODERM CGF DRESSING

APPLY A DRESSING THAT EXTENDS 3 CM BEYOND WOUND

Cleanse the wound site, rinse well and dry the surrounding skin. Note: As wound exudate is absorbed by the dressing, gel formation may be visible as a bubble on the outer surface of the dressing.

1. Remove the silicone release paper from the back of the dressing, minimizing finger contact with adhesive surface.
2. Gently roll dressing over wound; do not stretch.
3. Mold dressing into place with hand for secure adhesion around wound margins.
4. Secure edges of dressing with hypoallergenic tape if needed.

REMOVAL OF THE DRESSING

5. Press down on the skin with one hand and carefully lift an edge of the dressing with your other hand.

Maximum recommended wear time is up to seven days.

STORAGE PRECAUTIONS

Store at room temperature. Keep dry. Avoid refrigeration and exposure to high humidity.

This product is supplied sterile.

© 2021 ConvaTec Inc.

® indicates a trademark of ConvaTec Inc.

ConvaTec Inc.
Greensboro, NC 27409
USA
1-800-422-8811

ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

Sponsored in Australia by:
ConvaTec (Australia) Pty Ltd.
Unipark Monash, Building 2, Ground Floor
195 Wellington Road
Clayton, Victoria 3168
(P.O. Box 63, Mulgrave, Victoria, Australia 3170)

ConvaTec International Services GmbH
Mühlentalstrasse 36/38
8200 Schaffhausen, Switzerland

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