convatec

Aquacel™

Extra

Hydrofiber™ Dressing with Strengthening Fibre / Medicazione in Hydrofiber™ con Fibra Rinforzate / Pansement Hydrofiber™ avec Renfort Fibre / Apósito de Hydrofiber™ con fibra reforzante / Pensó Hydrofiber™ com Fibra de Reforço / Hydrofiber™ Verband mit Verstärkenden Fasern / Hydrofiber™ Förband med Förstärkande fibrer / Kultuvahvistettu Hydrofiber™-haavaside / Hydrofiber™ Bandage med förstärkade fibre / Hydrofiber™ Bandasje med Forsterkede Fiber / Ενδοαση Hydrofiber™ με ενωγυρτικές ίνες / مع تسليق موى Hydrofiber™ مع

INSTRUCTIONS FOR USE / ISTRUZIONI PER L'USO / NOTICE D'UTILISATION / INSTRUÇÕES DE USO / INSTRUÇÕES DE USO / GEBRUIKSIJNSTRUCTIES / BRUKSANVISNING / KÄYTTÖOHJEET / BRUGSANVISNING / BRUKSANVISNING / تعليمات الاستخدام /

ENGLISH

PRODUCT DESCRIPTION

Aquacel™ Extra™ Dressings are soft, sterile, dressings made from two layers of sodium carboxymethylcellulose stitched together with strengthening fibres. This dressing absorbs high amounts of wound fluid and bacteria and creates a soft, cohesive gel that intimately conforms to the wound surface, maintains a moist environment to support the body's healing process and aids in the removal of non-viable tissue from the wound (autolytic debridement).

INTENDED USE

Aquacel™ Extra™ Dressings have been designed to be used as a primary dressing under the direction of a healthcare professional as per the indications for use.

INTENDED USERS

Aquacel™ Extra™ Dressings are intended to be used by healthcare professionals, carers and patients under the direction of a Healthcare Professional.

TARGET PATIENT POPULATION

Aquacel™ Extra™ Dressings are intended to be used on patients with one of the wound types listed in the indications for use.

CLINICAL BENEFITS

Aquacel™ Extra™ Dressings absorb wound fluid and bacteria, providing a moist wound healing environment, aiding autolytic debridement and removing dead-space between the wound and dressing interface. Aquacel™ Extra™ Dressings are designed to manage excess exudate levels, which may further damage the wound bed and surrounding skin. Aquacel™ Extra™ Dressings provide a barrier to protect the wound bed from contamination.

INDICATIONS

Aquacel™ Extra™ Dressings are indicated for:

- Leg ulcers, including:
 - Venous ulcers
 - Arterial ulcers
- Diabetic ulcers including diabetic foot ulcers
- Pressure ulcers/injuries
- Surgical wounds
- Traumatic wounds
- Malignant wounds

CONTRAINDICATIONS

Aquacel™ Extra™ Dressings should not be used on individuals who are sensitive to or who had an allergic reaction to the dressing or any of its components.

PRECAUTIONS AND OBSERVATIONS

- This dressing is not intended for use as a surgical sponge.
- Sterility is guaranteed unless the pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or open prior to use and dispose of the device according to local regulations.
- Aquacel™ Extra™ Dressings are for single-use only and should not be re-used. Reuse may lead to increased risk of infection, cross-contamination and delayed healing.
- This dressing should not be used with other wound care products without first consulting a healthcare professional.
- Aquacel™ Extra™ Dressings are not compatible with petroleum-based products.
- Due to the sterilisation process, there may be a slight odor on opening the primary packaging.
- During the body's normal healing process, non-viable tissue is removed from the wound (autolytic debridement), which could initially make the wound appear larger.
- Newly formed blood vessels may occasionally produce bloodstained wound fluid following removal of the dressing.
- Aquacel™ Extra™ Dressings are designed to be MRI safe.
- A healthcare professional should be consulted if any of the following are observed during dressing changes, irritation (reddening, inflammation), maceration (whitening of skin), hyper-granulation (excess tissue formation), signs of infection (increased pain, bleeding, redness, tenderness of surrounding tissue, wound exudate), or a change in wound colour and/or odour.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local state and federal laws and regulations.
- If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the Competent Authority of the Member State in which the user and/or patient is established.

ASSOCIATED PRODUCT & THERAPIES

Aquacel™ Extra™ Dressings may be used alone or in combination with other wound care products such as; primary dressings, secondary dressings, securement bands, compression therapy devices or medicated topical treatments. Use of additional products and therapies should be undertaken with the direction of a Health Care Professional.

CONTACT AND DURATION OF USE

Aquacel™ Extra™ Dressings can be worn for up to 7 days, dressings should be changed earlier if clinically indicated.

DIRECTIONS FOR USE

- Before applying the dressing, cleanse the wound area with an appropriate wound cleanser.
- Aquacel™ Extra™ Dressings should overlap at least 1 cm onto the skin surrounding the wound.
- This dressing can be cut in any direction.
- When using Aquacel™ Extra™ Dressings in deep wounds, only fill the wound up to 80% as Aquacel™ Extra™ Dressings will expand to fill the wound space on contact with wound fluid.
- A secondary dressing is required to hold Aquacel™ Extra™ Dressings in place, such as:
 - a moisture retentive cover dressing such as DuoDERM™ Extra Thin in lighty to moderate exudate
 - a secondary cover dressing such as Aquacel™ Foam or ConvaMax™, in moderate to heavily exuding wounds;
 - for dry wounds refer to section FOR DRY WOUNDS below.
- See individual cover dressing package inserts for complete instruction for use.
- All wounds should be inspected regularly. Remove the Aquacel™ Extra™ Dressings when clinically indicated (i.e., leakage, excessive bleeding, increased pain) after a maximum of seven days.

FOR DRY WOUNDS

- Pre-moisten the Aquacel™ Extra™ Dressings by wetting with sterile saline prior to dressing the wound. The gelling properties of Aquacel™ Extra™ Dressings will help to maintain a moist wound and reduce the risk of maceration.
- Fully cover the dressing with a moisture retentive dressing such as DuoDERM™ Extra Thin to avoid drying out of the dressing and subsequent dressing adherence to the wound.
- If on removal, the dressing is dry, hydrate with sterile saline prior to removal to reduce risk of trauma.

Store in cool dry place. (10°C - 25°C/50°F – 77°F). Keep Dry.

If further information or guidance is needed, please contact Convatec Professional Services.

Made in UK

© 2025 Convatec

™/® indicates a trademark of the Convatec Group

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