**Hydroflex® Dressing with Silver Reinforced with Nylon / Hydroflex® Verband mit Silber verstärkt mit Nylon / Hydroflex® Vestir con Nylon reforzado con nylon / Pansament Hydroflex® à l'argent avec un renfort en Nylon / Apósito de Hydroflex® con plata reforzada con nylon / Hydroflex® vestir com Prata reforçado com Nylon / Hydroflex®verband mit silber, verstärkt mit nylon / Hydroflex® förband med silver, förstärkt med nylon / Hessa silberilläinen / Hydroflex®-pöytähoitoaine / Hydroflex® bandage med sølv, förstärkt med nylon / Hydroflex® nylonforstärkt bandage med sølv / Apósito de plata reforzada con nylon / Hydroflex® Apósito en nylon / Hydroflex® مع مادة النايون / Hydroflex® مع مادة النايون / Nylon Taktiylit, Gümiyle İyelen Hydroflex®® Tara Otası**

**INSTRUCTIONS FOR USE / GEBRUCHSANWEISUNG / ISTRUZIONI PER L'USO / NOTICE D'UTILISATION / INSTRUCCIONES DE USO / INSTRUCÇÕES DE USO / BRUKSANVISNING / BRUKSANVISNING / KÄTTÖOHJEET / BRUGSVEJLEDNING / BRUKSANVISNING / OAHITSE XPIRE / تاریخ استعمال / KULLANMIN TALIMALLAR**

**PRODUCT DESCRIPTION**  
**AQUACEL™ Ag BURN Hydroflex®** Dressing with Silver Reinforced with Nylon, silver-impregnated, antimicrobial, absorbent dressing is a soft, sterile, non-woven hydroentangled dressing comprised of Hydroflex® Technology (sodium carboxymethylcellulose) with nylon thread used to stitch bond the product along the length. The dressing contains 1.2% w/w ionic silver. The silver in the dressing kills a broad spectrum of wound bacteria, and aids in creating an antimicrobial environment. This dressing absorbs high amounts of wound fluid and bacteria and creates a soft, cohesive gel that intimately conforms to the wound surface, maintains a moist environment and aids in the removal of non-viable tissue from the wound (autolytic debridement). A moist wound environment and control of wound bacteria supports the body's healing process and helps reduce the risk of wound infection.

**INTENDED USE**  
AQUACEL™ Ag BURN dressing is intended to be used in the management of burns.

**INTENDED USERS**  
AQUACEL™ Ag BURN is intended to be used by healthcare professionals.

**TARGET PATIENT POPULATION**  
Patients with first- and second-degree burns.

**CLINICAL BENEFIT**  
The intended clinical benefit of AQUACEL™ Ag BURN is to provide wound healing, infection management, pain reduction and improved wear time.

**INDICATIONS**  
AQUACEL™ Ag Burn may be used for the management of superficial (first degree) and partial thickness (second degree) burns.

**CONTRAINDICATIONS**  
Do not use on full thickness (third degree) burns. The dressing should not be used on individuals who are sensitive to or have a known allergy to silver, sodium carboxymethylcellulose or nylon.

**PRECAUTIONS AND OBSERVATIONS**  
• Do not reuse. These devices are for single use only and should not be re-used. Reuse may lead to increased risk of infection and cross-contamination and delayed healing.  
• Dressings should not be used if the pouch is damaged or opened prior to use as sterility cannot be guaranteed. Due to the sterilization process there may be slight odour on opening the primary packaging.  
• The dressing should be inspected at each dressing change and appropriate treatment given if there are signs of infection in the wound/peri-wound area, or if there is an adverse reaction to the dressings.

• In the management of partial thickness hand burns, adherence to the wound bed of the AQUACEL™ Ag BURN dressing, is a desired characteristic.  
• If the dressing is difficult to remove, it should be soaked with sterile water or sterile saline until it can be easily removed.

• Because the AQUACEL™ Ag BURN dressing provides a moist environment to support the growth of new blood vessels, the delicate newly formed blood vessels may be occasionally produce blood stained wound fluid.  
• AQUACEL™ Ag BURN dressing is not compatible with oil-based products, such as petroleum.

• AQUACEL™ Ag BURN dressing is not intended for use as a surgical sponge. Silver-containing products may cause skin discoloration following prolonged use; normal use of the dressing can continue unless wound healing is affected.

Clinicians/Healthcare professionals should be aware that there are very limited data on prolonged and repeated use of silver containing products, particularly in children and neonates.

AQUACEL™ Ag BURN Dressings do not replace the need for systemic or other adequate infection treatment.

**DIRECTIONS FOR USE**  
For partial thickness burns, the AQUACEL™ Ag BURN Dressings may be left in place for up to 14 days before changing unless an earlier change is clinically indicated. If the wound is infected, frequent inspection of the wound is recommended. Silvered dressings, such as the AQUACEL™ Ag BURN Dressing,

are intended for short-term use and as such consecutive use of Aquacel™ Ag BURN Dressings should not exceed 28 days.  
**Applying the dressing**  
• Before applying the dressing, cleanse the wound area with an appropriate wound cleanser.  
• The AQUACEL™ Ag BURN Dressing should overlap at least 2cm onto the skin surrounding the burn and any other adjacent AQUACEL™ Ag BURN Dressing.

• The AQUACEL™ Ag BURN Dressing may be secured, with a secondary dressing or retention bandage. If a secondary dressing is used, remove regularly to inspect the AQUACEL™ Ag BURN Dressing, which can remain in place on the wound.  
• All dressings should be inspected frequently.  
• Dressings should be changed as clinically indicated.  
• For partial thickness burns AQUACEL™ Ag BURN Dressing may be left in place for up to 14 days or until clinically indicated. If the wound is infected, frequent inspection of the wound may be necessary.

**ASSESSMENT OF ADHERENCE**  
The AQUACEL™ Ag BURN Dressing should be inspected for adherence to the burn. Remove the secondary cover dressing and inspect the dressing for possible adherence.  
• If the dressing is saturated or leaking it should be easily removed and replaced, otherwise it can be left in place.  
• If the dressing is adhered and there appears to be no signs of infection the AQUACEL™ Ag BURN dressing can be left in place for up to 14 days and will act as a synthetic scab, debriding as a natural scab would.

**Removal of the dressing**  
As the wound re-epithelializes, the AQUACEL™ Ag BURN Dressing should detach or be easily removed. If you have difficulty removing the dressing, it should be soaked with sterile water or sterile saline until it can be easily removed.

Store at room temperature (10°C- 25°/50°-77°F). Keep dry.  
After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local laws and regulations.

If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national competent authority.

https://ec.europa.eu/tools/euadmed. This is the Summary of Safety and Clinical Performance (SSCP) location after the launch of European Database on Medical Devices (EUADMED).

If further information or guidance is needed, please contact ConvaTec Professional Services.  
© 2022 ConvaTec Inc.  
\*\*/™ indicates a trademark of ConvaTec Inc.

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\*\*/™ indique une marque déposée de ConvaTec Inc.