<p>Sterile / Steriel / Steriili / Steril</p>

<p>A Tapered Hydrocolloid Dressing with a Change Indicator / Hydrokolloidverband mit Verbandwechsellindikator / MedicaZIONE idrocolloidaLe con Indicatore per il Cambio / Pansement hydrocolloïde à bords effilés avec indicateur de changement / Een hydrocolloid verband met aflopende randen en verwisselindicator / Ett Hydrokolloidförband med avfasade kanter och bytesindikator / Reunolita ohennettu hydrokolloidisidos, joka näyttäää, milloin sidos tulee vaihtaa / En hydrokolloid förbinding med affaset kant og skifteindikator / Hydrokolloid sårbandasje med avfasede kanter og skifteindikator / 帶有纖薄黏貼邊緣及更換指示的水膠體敷料</p>
<p>INSTRUCTIONS FOR USE / GEBRAUCHSANWEISUNG / ISTRUZIONI PER L'USO / MODE D'EMPLOI / GEBRUIKSAANWIJZING / BRUKSANVISNING / KÄYTTÖOHJEET / BRUGSANVISNING / BRUKSANVISNING / 使用方法</p>

ENGLISH

PRODUCT DESCRIPTION

DuoDerm Signal™ is an adhesive, bordered hydrocolloid wound dressing with a change indicator and tapered edges. An indicator (a printed green line which marks the outer limit for signs of exudate spread prior to dressing change) helps to determine when to change the dressing. The wound fluid appears as a bubble moving towards the indicator line. The dressing is tapered to be thinner at the edges than at the centre, to reduce the incidence of rolling up during use. The adhesive layer forms a cohesive gel when in contact with wound exudate.

The dressing absorbs wound fluid and creates a moist environment which supports the body's healing process and aids in the removal of non-viable tissue from the wound (autolytic debridement), without damaging newly formed tissue.

DuoDerm Signal may be used as a primary dressing or as a secondary dressing. It may be used alone or in combination with other wound care products.

DuoDerm Signal acts as a barrier to the wound against bacterial, viral and other external contamination. It minimizes the potential of exposure to nosocomial or infectious agents. The dressing has been shown in laboratory experiments (7 day testing period) to block the passage of bacteria and the Phi-X 174 virus (a validated surrogate for the Human Immunodeficiency Virus (HIV-1) and Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV)) provided the dressing remains intact and there is no leakage. The use of this dressing neither guarantees nor warrants against AIDS, HBV or HBC transmission.

INDICATIONS

Under the supervision of a healthcare professional, DuoDerm Signal may be used for wounds such as:

- leg ulcers (venous stasis ulcers, arterial ulcers and leg ulcers of mixed aeti-ology), diabetic ulcers and pressure ulcers/sores (partial & full thickness)
- surgical wounds (post-operative lefts to heal by secondary intention, donor sites, dermatological excisions)
- second degree burns
- traumatic wounds including:
 - abrasions
 - lacerations
 - minor cuts

KONTRAINDICATIONS

DuoDerm Signal should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.

PRECAUTIONS AND OBSERVATIONS

- Sterility is guaranteed unless immediate product packaging is damaged or opened prior to use. This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use.
- Appropriate supportive measures should be taken where indicated (e.g. use of graduated compression in the management of ve- nous leg ulcers or pressure relief measures in the management of pressure ulcers).
- The control of blood glucose, as well as appropriate supportive measures, e.g. pressure off-loading, should be provided with diabetic foot ulcers.
- During the body's normal healing process, non-viable tissue (autolytic de-bridement) is removed from the wound, which could make the wound ap-pear larger after the first few dressing changes. If the wound continues to get larger after the first few dressing changes, consult a healthcare profes- sional.
- Should you observe irritation (reddening, inflammation), maceration (whitening of skin), hypergranulation (excess tissue formation) or sensitivity (allergic reaction) consult a healthcare professional.
- Frequent dressing changes on wounds with damaged or delicate skin surrounding the wound is not recommended.
- Colonization of chronic wounds is common and is not a contraindication to the use of the dressing. Further, the dressing may be used on infected wounds under medical supervision together with appropriate therapy and frequent monitoring of the wound.
- The use of this dressing neither guarantees nor warrants against AIDS, HBV or HCV transmission.

DIRECTIONS FOR USE

1. Wound Site Preparation and Cleansing

Before using the dressing, clean the wound with an appropriate wound cleansing agent or saline and dry the surrounding skin.

2. Dressing Preparation and Application

- Choose a dressing size and shape to ensure that the dressing is 3cm larger than the wound area in any one direction.
- Remove the release paper from the back being careful to minimize finger contact with the adhesive surface.
- Hold the dressing over the wound and line up the centre of the dress- ing with the centre of the wound. Place the dressing directly over the wound.
- For difficult to dress anatomical locations, such as heels or the sacrum, please use the appropriate shape dressing.
- Discard any unused portion of the product after dressing the wound.

3. Dressing Change and Removal

- NOTE:** The dressing should be inspected frequently (depending on the level of exudate) to see whether any part of the bubble has reached the change indicator, or whether there has been any leakage. If any of these occur, the dressing should be changed as clinically indicated. As wound fluid is absorbed by the dressing, gel formation may be visible through the dressing.
- The dressing should be changed when: clinically indicated, the bubble has reached the indictor or when used up to a maximum of seven days.
- Press down gently on the skin and carefully lift one corner of the dress- ing until it is no longer adhered to the skin. Continue until all edges are free. Carefully lift away the dressing.

Store at room temperature (10° C - 25° C/50-77° F). Keep pack dry when storing.

If further information or guidance is needed, please contact ConvaTec Professional Services.

© 2021 ConvaTec Inc.
*/™ indicates a trademark of ConvaTec Inc.

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*/™ indikerer varemærker som tilhører ConvaTec Inc.