CarboFlex®

Odour Control Dressing / Geruchsneutralisierende Wundauflage / Medicazione per il controllo degli odori / Pansement destiné au contrôle des odeurs / Apósito de control del olor / Penso para Controlo do Odor / Geurbestrijdingsverband / Luktabsorberande förband / Hajua imevã sidos / Lugtabsorberende bandage / Luktabsorberende bandasje / Επίθεμα Ελέγξου Οσμών / منسدة للتحكم بالرائحة

en ENGLISH

PRODUCT DESCRIPTION
CarboFlex® Odour Control Dressing is a sterile five-layer non-adhesive dressing with an absorbent wound contact layer (containing Alginate and Hydrofibre™) an activated charcoal central pad, an absorbent padding layer and soft flocking outer contained within a peel pouch.

INTENDED USE
CarboFlex® Odour Control Dressing is intended for the management of malodorous wounds. It may be used as a primary or secondary dressing as per the indications for use.

INTENDED USERS
CarboFlex® Odour Control Dressing is intended to be used by healthcare professionals and by carers and patients under the direction of a Healthcare Professional.

TARGET PATIENT POPULATION
CarboFlex® Odour Control Dressing is intended for use in adults and children. It is not recommended for infants less than 12 months of age.
CarboFlex® is intended to be used on patients with one of the wound types listed in the indications for use.

CLINICAL BENEFITS
The clinical benefits of CarboFlex® Odour Control Dressing are:

- Management of malodorous wounds
- Absorb and control exudate
- To support maintaining a moist environment for wound healing

INDICATIONS
CarboFlex® Odour Control Dressing is indicated for the management of pressure ulcers, venous leg ulcers, arterial leg ulcers, mixed aetiology ulcers and diabetic foot ulcers.

CONTRAINDICATIONS

CarboFlex® Odour Control Dressing should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or any of its components including calcium sodium alginate or other seaweed derivatives.

PRECAUTIONS AND OBSERVATIONS

- The use of CarboFlex® Odour Control Dressings is not recommended for infants less than 12 months of age due to the presence of manganese within the dressing (≤ 20ppm). Rarely, irritation and/or maceration of the skin and hypergranulation have been reported.
- Sterility is guaranteed unless the pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or open prior to use and dispose of the device according to local regulations.
- Do not use after stated expiry date.
- The dressing may be used on infected wounds only under the care of a healthcare professional.
- CarboFlex® Odour Control Dressing is not compatible with petroleum-based products.
- CarboFlex® Odour Control Dressing is not designed to be cut.
- CarboFlex® Odour Control Dressing is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use.
- Should you observe irritation (reddening, inflammation), maceration (whitening of skin), hypergranulation (excess tissue formation) or sensitivity (allergic reaction), consult a healthcare professional.
- After use, this product may be a potential biohazard. Discard as per local clinical protocols.
- If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

ASSOCIATED PRODUCTS & THERAPIES

CarboFlex® Odour Control Dressing may be used as a primary or secondary dressing. It may be used alone or in combination with other wound care products and therapies such as but not limited to: compression therapy, primary dressings and wound fillers as directed by the healthcare professional.

CONTACT & DURATION OF USE
CarboFlex® Odour Control Dressing can be worn for up to 3 days. Dressings should be changed earlier if clinically indicated. It is safe to continue using CarboFlex® Odour Control Dressing for as long as clinically indicated.

DIRECTIONS FOR USE

Wound site preparation and skin cleansing:

These dressings are sterile and should be handled appropriately. If required, the wound should be debrided and necrotic tissue removed. Cleanse the wound site, rinse well and dry the surrounding skin.

Dressing preparation and application:

Do not cut dressing.

- Choose a dressing sized so that it is larger than the wound edge leaving a sufficient border. For shallow wounds the dressing may be placed directly onto the wound as a primary dressing and for cavity wounds CarboFlex® can be laid over a wound filler/gel as a secondary dressing.
- Place fibrous (non-shiny) surface on wound with primary dressing.
- Secure CarboFlex® in place over the wound with tape or other appropriate material. The absorbent wound contact layer will take up exudate and form a soft gel.
- Change the dressing when clinically indicated or if there is exudate strike through of the top layer or if the odour is no longer being absorbed. With noninfected malodorous wounds, CarboFlex® may be left undisturbed for up to three days.

Removal and disposal

- Carefully lift away from the wound by grasping one corner of the dressing. Removal from lightly exuding wounds may be assisted by saturating the dressing with sterile normal saline.
- After use, this product may be a potential biohazard. Discard as per local clinical protocols.

STORAGE PRECAUTIONS

Store at room temperature (10-25°C/50-77°F). Avoid refrigeration and exposure to high humidity. Keep Dry. If further information or guidance is needed, please contact Convatec Professional Services.

© 2025 Convatec

™/® indicates a trademark of the Convatec Group

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