
ESENTA™

Stomahesive®

Paste / Pasta / Pasta

Protective Skin Barrier & Filler / Huidbeschermers & Opvulmiddel / Koruyucu Cilt Bariyeri ve Dolgusu

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en ENGLISH

Intended Use:

Protect exposed skin between the base of the stoma and the opening in the wafer, and as a filler for skin folds and uneven skin surfaces.

Indications for Use: For use with a colostomy, ileostomy or urostomy.

PRODUCT DESCRIPTION

Stomahesive Paste is a convenient, highly effective filler, sealant, and protective skin barrier specifically designed for use around colostomy, ileostomy, or urinary stomas.

CONTRAINDICATION

Do not use on individuals with a known sensitivity to the paste or its ingredients.

IMPORTANT:

1. For external use only. Do not apply into an open wound or skin ulcer.
2. Contains alcohol.
3. Should irritation occur, discontinue use and consult a health care professional.

DIRECTIONS FOR USE

APPLICATION TO SKIN BARRIER

- Prepare the Flange/Skin Barrier as described in the Instructions for Use accompanying the Flange/Skin Barrier.
- Apply a ring of paste on the back of the Flange/Skin Barrier around the hole. Do not spread the paste. Wait one minute before applying the Flange/Skin Barrier.
- For best results, squeeze tube from bottom and flatten as you go up.

APPLICATION TO SKIN

Apply Stomahesive Paste to fill depressions and irregular skin contours. If required, use a moistened finger or spatula to spread the paste evenly. (1)

GENERAL INSTRUCTIONS

- Apply the Flange/Skin Barrier and press down firmly

around the base of the stoma. (2)

- Additional paste may be applied to further protect the skin around the stoma at any time.
- Tighten cap securely.

REMOVING STOMAHESIVE PASTE

Stomahesive Paste is easily removed after it has been in place for 48 hours. If the paste must be removed sooner, use warm water. It is not necessary to remove all traces of the paste prior to a new application.

DISPOSAL

If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

Handle and dispose of in accordance with accepted medical practice and applicable local laws and regulations.

For more information on Convatec products, contact us at the address below.

© 2025 Convatec Inc.

®/™ indicates a trademark of Convatec Inc.

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®/™ Convatec Inc. firmasının tescilli ticari markasıdır.
Convatec Limited
UNOMEDICAL A/S