FoamLite™

Adhesive Adhésif

Adhesive Silicone Foam Dressing
Aposito de Espuma con Adhesivo de Silicona
Pansement hydrocellulaire adhesif silicone
Curativo Adesivo de Espuma e Silicone

en ENGLISH

Product Description

Foam Lite™ Convatec dressing is a thin, conformable, laminated polyurethane layered foam dressing. The dressing consists of a waterproof outer polyurethane film and an absorbent foam pad with a perforated gentle silicone adhesive.

The outer film layer provides a waterproof, viral and bacterial barrier which protects the wound from external contamination helping reduce the risk of infection. The dressing film backing acts as a barrier to the wound against bacteria and blood borne viral pathogens.

The wound contact surface is coated with a perforated gentle silicone adhesive layer which helps provide secure skin friendly adhesion.

Foam Lite™ Convatec dressing is designed for low exuding wounds and helps maintain a moist wound environment which can help promote faster wound closure. Foam Lite™ Convatec dressing is designed to remain in place without the need for secondary retention. Foam Lite™ Convatec dressing can be repositioned during application.

Indications

Foam Lite™ Convatec dressings are intended for the management of a wide range of low / non exuding wounds. Under medical supervision Foam Lite™ Convatec dressings may be used for the management of both chronic and acute wounds such as, but not limited to:

- Leg ulcers, pressure ulcers and diabetic ulcers;
- Surgical wounds (e.g. post operative wounds left to heal by secondary intent and donor sites);
- Partial thickness (second degree) burns;
- Traumatic wounds (e.g. abrasions, lacerations, blisters, minor cuts and skin tears).

Contraindications

Foam Lite™ Convatec dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.

Foam Lite™ Convatec dressings are not compatible with oil based products or emollients such as petrolatum.

Do not use Foam Lite™ Convatec dressings with oxidising agents such as hydrogen peroxide or hypochlorite solutions.

Precautions and Observations

Caution: Sterility is guaranteed unless pouch is damaged or opened prior to use. This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use. This wound dressing should not be used with other wound care products without first consulting a healthcare professional.

During the body's normal healing process, non-viable tissue may be removed (autolytic debridement) from the wound, which could make the wound appear larger after the first few dressing changes. Should you observe irritation (reddening, inflammation), maceration (whitening of skin), or hypergranulation (excess tissue formation), consult a healthcare professional. The wound should be inspected during dressing changes for: (1) signs of infection (increased pain, bleeding, warmth/redness of surrounding tissue, wound exudate), (2) there is a change in wound color and/or odor, (3) any other unexpected symptoms occur (e.g. maceration or hypergranulation).

The use of Foam Lite™ Convatec dressings have not been studied in wounds due to herpes simplex or impetigo. In addition, for leg ulcers, pressure ulcers, diabetic ulcers, second degree burns, donor sites, surgical or traumatic wounds left to heal by primary or secondary intention:

- Treatment of the wound types listed above should only be under the supervision of a healthcare professional.
- Appropriate supportive measures should be taken where indicated (e.g. use of graduated compression bandaging in the management of venous leg ulcers or pressure relief measures in the management of pressure ulcers).
- Colonization of chronic wounds is common and is not a contraindication to the use of the dressing. The dressing may be used on infected wounds under medical supervision together with appropriate therapy and frequent monitoring of the wound.
- The control of blood glucose, as well as appropriate supportive measures, should be provided with diabetic foot ulcers.
- In second degree burns, consider alternate surgical procedures if the wound has not reepithelialized after 14 days.

Foam Lite™ Convatec dressings may be cut to shape and size as directed by a healthcare professional. If the dressing is cut, additional tape may be required to secure it in place. If the dressing is cut, the viral / bacterial barrier provided by the dressing may be compromised.

Directions for Use

1. If the immediate sterile pouch is damaged do not use.

2. Wound Site Preparation and Cleaning:

Before applying the dressing, cleanse the wound area with appropriate wound cleanser in accordance with normal procedure and dry the surrounding skin. Ensure that the surrounding skin is free from any emollients or oil based products that might affect dressing adhesion.

3. Dressing Preparation and Application

- Select the appropriate dressing size – choose a dressing size and shape to ensure that the central pad is at least 1cm (1/3 inch) larger than the wound area.
- Remove the dressing from the sterile pack, minimising finger contact with the wound contact surface of the dressing and remove the release liner.
- Hold the dressing over the wound and line up the centre of the dressing with the centre of the wound. Smooth the dressing over the wound and ensure that the dressing is adhered all around the wound.
- Foam Lite™ Convatec dressings can be cut especially to dress wounds on the heel, elbows and other awkward areas. Additional tape or other means of fixation will be required to secure the dressing in place.
- Discard any unused portion of the product after dressing the wound.

4. Dressing Removal

- The dressing should be changed when clinically indicated (e.g. leakage, bleeding, increased pain, suspicion of infection). Maximum recommended wear time is up to seven days.
- To remove the dressing, press down gently on the skin and carefully lift one corner of the dressing. Continue until all edges are free. Carefully lift away the dressing and discard according to local clinical protocols.

Store at room temperature (10° C - 25° C/50° F - 77° F).
Protect from light. Keep dry.

If further information or guidance is needed, please contact Convatec Professional Services or visit us online at www.convatec.com.

Made in UK

© 2026 Convatec

™/® indicates a trademark of the Convatec Group

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Convatec (Australia) Pty Ltd.
Getz Bros. Philippines, Inc.
Importado por: Convatec PERU S.A.C.