INSTRUCTIONS FOR USE / GEBRAUCHSANLEITUNG / ISTRUZIONI PER L'USO / CONSEILS D'UTILISATION / INSTRUCCIONES DE USO / INSTRUCÇÕES DE USO / GEBRUIKSIJNSTRUCTIE / BRUKSANVISNING / RÄTTORDBUEK / BRUGSANVISNING / BRUKSANVISNING / LEIBENINGANGUR UN NOTKVAÐAÐÆT XPHÆSE / تعليمات الإستعمال

en ENGLISH

PRODUCT DESCRIPTION
AQUACEL™ Ribbon Dressing with Strengthening Fibre is a soft, sterile, non-woven ribbon dressing composed of sodium carboxymethylcellulose and regenerated cellulose fibre for strengthening. This comfortable and highly absorbent dressing absorbs wound fluid and transforms into a soft gel, which maintains a moist environment to support the body's healing process and aid in the removal of nonviable tissue from the wound (autolytic debridement), without damaging newly formed tissue.

INTENDED USE

AQUACEL™ Ribbon Dressings have been designed to be used as a primary dressing for cavity wounds. They are intended to be used under the direction of a healthcare professional as per the indications for use.

INTENDED USERS
AQUACEL™ Ribbon Dressings are intended to be used by healthcare professionals, carers and patients under the direction of a Healthcare Professional.

TARGET PATIENT POPULATION
AQUACEL™ Ribbon Dressings are intended to be used on patients with one of the wound types listed in the indications for use.

CLINICAL BENEFITS

AQUACEL™ Ribbon Dressings absorb wound fluid and bacteria, providing a moist wound healing environment, aiding autolytic debridement and removing dead space between the wound and dressing surface. AQUACEL™ Ribbon Dressings are designed to manage exudate levels which may further damage the wound bed and surrounding skin. AQUACEL™ Ribbon Dressings provide a barrier to protect the wound bed from contamination.

INDICATIONS

AQUACEL™ Ribbon Dressings are indicated for:

CONTRAINDICATIONS
AQUACEL™ Ribbon Dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or any of its components.

PRECAUTIONS AND WARNINGS

• This dressing is not intended for use as a surgical sponge.
• Sterility is guaranteed unless pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or open prior to use and dispose of the device accordingly to local regulations.
• AQUACEL™ Ribbon Dressings are for single-use only and should not be re-used. Re-use may lead to increased risk of infection, cross-contamination or delayed healing.
• This dressing should not be used with other wound care products without first consulting a healthcare professional.

• AQUACEL™ Ribbon Dressings can be cut, if required prior to use, across the narrower dimension of the ribbon. Cutting along the longer dimension of the ribbon should not be performed.

• AQUACEL™ Ribbon Dressings are not compatible with petroleum-based products.
• Due to the sterilisation process there may be a slight odour on opening the primary packaging.
• During the body's normal healing process, non-viable tissue is removed from the wound (autolytic debridement), which could initially make the wound appear larger. As the wound forms blood vessels, it occasionally produce bloodstained wound fluid following removal of the dressing.

• AQUACEL™ Ribbon Dressings are designed to be MRI safe.
• A healthcare professional should be consulted if any of the following are observed during dressing changes, irritation (redness, inflammation), maceration (whitening of skin), hypergranulation (excess tissue formation), signs of infection (increased pain, bleeding, warmth/redness of surrounding tissue, wound exudate), or a change in wound colour and/or odour.
• After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local state and federal laws and regulations.
• During the use of this device as a result of its use a serious incident has occurred, please report it to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

ASSOCIATED PRODUCTS & THERAPIES

AQUACEL™ Ribbon Dressings may be used alone or in combination with other wound care products such as primary dressings, secondary dressing, securement bandages, compression therapy devices or medicated topical treatments. Use of additional products and therapies should be undertaken with the supervision of a Health Care Professional.

CONTACT AND DURATION OF USE
AQUACEL™ Ribbon Dressings can be worn for up to 7 days, dressings should be changed either if clinically indicated.

DIRECTIONS FOR USE

• Before applying the dressing, cleanse the wound area with an appropriate wound cleanser.
• AQUACEL™ Ribbon Dressing with Strengthening Fibre should overlap at least 1 cm onto the skin surrounding the wound.
• When using AQUACEL™ Ribbon Dressing with Strengthening Fibre in deep cavity wounds, leave at least 2.5 cm outside the wound for easy retrieval.
• Only pack deep wounds up to 80%, as AQUACEL™ Ribbon, Dressing with Strengthening Fibre will expand to fill the wound, spaces on contact with wound fluid.
• Apply the dressing to the wound with a moisture retentive dressing such as DuoDERM™ Extra Thin™, or other appropriate secondary dressing. See relevant instruction package inserts for complete instructions for use.
• All wounds should be inspected frequently.
• Remove the AQUACEL™ Ribbon Dressing with Strengthening Fibre when clinically indicated (i.e., leakage, excessive bleeding, or suspicion of infection).
• AQUACEL™ Ribbon Dressing with Strengthening Fibre can be left in place for up to 7 days, where clinically indicated.


FOR DRY WOUNDS

In addition to the directions for use set forth above:
• Place the AQUACEL™ Ribbon Dressing with Strengthening Fibre on the wound and wet with sterile water or saline over the wound area only.
• To cover the dressing with a moisture retentive dressing such as DuoDERM™ Extra Thin to avoid drying out of the dressing and subsequent dressing adherence to the wound.

Discard any unused portion of the dressing
• If the dressing or packaging is damaged, do not use.
• Store in a cool dry place.
• If there is information or guidance is needed, please contact Convatec Professional Services.

© 2023 Convatec Inc.

*™ indicates a trademark of Convatec Inc.

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