R ONLY  
STERILE/EO

**convatec**

# Aquacel™ Ag

Foam

**INSTRUCTIONS FOR USE / INSTRUCCIONES DE USO / MODE D'EMPLOI / INSTRUÇÕES DE USO / 自黏性及非自黏性含銀親水性纖維泡棉敷料 / 使用说明书 / 사용설명서**

**Aquacel™ Ag foam adhesive and non adhesive Hydrofiber™ foam dressing with silver / Aquacel™ Ag foam adesivo y no adhesivo Apósito de foam de Hydrofiber™ con plata / Pansement mousse Aquacel™ Ag adhésif et non adhésif Pansement mousse Hydrofiber™ avec argent / Aquacel™ Ag foam adesivo e não adesivo Curativo de espuma de Hydrofiber™ com prata / 術部器輪字第027127號 使用前請務必詳閱原廠之使用說明書並遵照指示使用。 / 羧甲基纤维素钠银泡沫敷料 (粘性和非粘性) / 접착성 & 비접착성Aquacel™ Ag 폼 은 함유 Hydrofiber™ 폼 드레싱**

EN **ENGLISH**

#### PRODUCT DESCRIPTION

Aquacel™ Ag foam dressing adhesive and non adhesive, is a range of sterile Hydrofiber™ foam wound dressings consisting of a waterproof outer polyurethane film and a multi-layered absorbent pad, with the adhesive dressing having a silicone adhesive border. The multi-layered absorbent pad contains a layer of polyurethane foam and a non-woven wound contact layer of Hydrofiber™ technology (sodium carboxymethylcellulose). The dressing's Hydrofiber™ technology wound contact layer contains 1.2 % w/w ionic silver.

The silver in the dressing kills a broad spectrum of wound bacteria held in the dressing. This dressing absorbs high amounts of wound fluid and bacteria and creates a soft, cohesive gel that intimately conforms to the wound surface, maintains a moist environment and aids in the removal of non-viable tissue from the wound (autolytic debridement) without damaging healthy tissue. The moist environment is known to support the body's healing process. The outer film layer provides a waterproof, viral and bacterial barrier which protects the wound from external contaminants. The film also helps to manage the moisture vapor transmission of the exudate absorbed by the dressing.

The adhesive dressing range has a silicone border which provides secure, skin-friendly adhesion.

Aquacel™ Ag foam dressing may be used as a primary or as a secondary dressing. It may be used alone or in combination with other wound care products and can be cut to shape and size, as directed by a healthcare professional. If the adhesive dressing is cut, additional tape may be required to secure it in place. Additional tape, or other means of fixation, will be required to secure the Aquacel™ Ag foam non adhesive dressing in place.

Aquacel™ Ag foam adhesive dressing film backing acts as a barrier to bacteria and viruses.

Aquacel™ Ag foam dressing contains silver as a preservative to prevent or minimize microbial growth within the dressing.

#### INDICATIONS

Aquacel™ Ag foam dressings may be used for the management of both chronic and acute wounds, such as:

- Partial thickness (second degree) burns;
- Diabetic foot ulcers, leg ulcers, (venous stasis ulcers, arterial ulcers and leg ulcers of mixed etiology) and pressure ulcers/sores (partial & full thickness);
- Surgical wounds;
- Traumatic wounds;
- Wounds that are prone to bleeding, such as wounds that have been mechanically or surgically debrided and donor sites;
- Abrasions;
- Lacerations;
- Minor cuts;
- Minor scalds and burns.

#### CONTRAINDICATIONS

Aquacel™ Ag foam dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.

#### PRECAUTIONS AND OBSERVATIONS

**Caution:** Sterility is guaranteed unless pouch is damaged or opened prior to use.

This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use.

This wound dressing should not be used with other wound care products without first consulting a healthcare professional.

During the body's normal healing process, non-viable tissue may be removed (autolytic debridement) from the wound, which could make the wound appear larger after the first few dressing changes.

Should you observe irritation (reddening, inflammation), maceration (whitening of skin), or hypergranulation (excess tissue formation), consult a healthcare professional.

The wound should be inspected during dressing changes for (1) signs of infection (increased pain, bleeding, warmth/redness of surrounding tissue, wound exudate), (2) there is a change in wound color and/or odor, (3) any other unexpected symptoms occur (e.g. maceration or hypergranulation).

The dressing may be used on infected wounds only under the care of a healthcare professional.

**The use of Aquacel™ Ag foam dressing adhesive and non adhesive has not been studied in wounds due to herpes simplex or impetigo.**

While it has been observed that certain silver-containing products may cause transient skin discoloration following prolonged use, clinical studies with products containing Hydrofer™ technology with ionic silver have shown no such skin discoloration.

**In addition, for leg ulcers, (venous stasis ulcers, arterial ulcers and leg ulcers of mixed etiology), diabetic ulcers, pressure ulcers/sores, partial thickness (second degree burns), donor sites and surgical, traumatic wounds left to heal by primary or secondary intention:**

- Treatment of the wound types listed above should only be under the supervision of a healthcare professional.
- Appropriate supportive measures should be taken where indicated (e.g. use of graduated compression bandaging in the management of venous leg ulcers or pressure-relief measures in the management of pressure ulcers/sores, systemic antibiotics and frequent monitoring in the treatment of wound infection, control of blood glucose for diabetic ulcers, etc.).
- In partial thickness (second degree) burns, consider alternate (surgical) procedures if the wound has not been reepithelialized after 14 days.
- Colonization of chronic wounds is common and is not a contraindication to the use of the dressing.
- Aquacel™ Ag foam dressing adhesive and non adhesive was determined to be MR-Safe according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. *Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.*

#### DIRECTIONS FOR USE

**1. If the immediate sterile product pouch is damaged, do not use.**

**2. Wound Site Preparation and Cleaning:** Before applying the dressing, cleanse the wound area with an appropriate wound cleanser and dry the surrounding skin.

**3. Dressing Preparation and Application:**

- Choose a dressing size and shape to ensure that the central absorbent pad (area within the adhesive window), is 1cm (1/3 inch) larger than the wound area.
- Remove the dressing from the sterile pack, minimizing finger contact with the wound contact surface and the adhesive surface where applicable. Remove the release liner if using the adhesive dressing.
- The dressing can be cut to shape for convenience.
- Hold the dressing over the wound and line up the center of the dressing with the center of the wound. Place the pad directly over the wound. For the adhesive dressing, smooth down the adhesive border.
- An appropriate retention bandage or tape should be used to secure the dressing in place if the dressing does not have an adhesive border or if the adhesive dressing has been cut.
- For difficult to dress anatomical locations, such as the heel or the sacrum, the specially shaped adhesive dressings may be used.
- Discard any unused portion of the product after dressing the wound.

**4. Dressing Removal:**

- The dressing should be changed when clinically indicated (i.e.: leakage, bleeding, increased pain, suspicion of infection). Maximum recommended wear time is up to seven days.
- The wound should be cleansed at appropriate intervals.
- To remove the dressing, press down gently on the skin and carefully lift one corner of the dressing. Continue until all edges are free. Carefully lift away the dressing and discard according to local clinical protocols.

Store at room temperature (10° C- 25° C/50° F- 77° F). Protect from light. Keep dry.

If further information or guidance is needed, please contact ConvaTec Professional Services or visit us online at [www.convatec.com](#).

Made in UK

© 2022 ConvaTec Inc.
™/™ indicates a trademark of ConvaTec Inc.
Aquacel, the Aquacel logo and Hydrofiber are registered trademarks in the United States.

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UNOMEDICAL A/S
ConvaTec Limited,
ConvaTec (Australia) Pty Ltd.,