
EsteemBody™

(en) INSTRUCTION FOR USE	(da) BRUGSANVISNING	(sr) UPUTSTVO ZA UPOTREBU
(de) GEBRAUCHSANWEISUNG	(no) BRUKSANVISNING	(lv) LIETOŠANAS PAMĀCĪBA
(it) ISTRUZIONI PER L'USO	(el) Οδηγίες χρήσης	(lt) NAUDOJIMO INSTRUKCIJA
(fr) NOTICE D'UTILISATION	(hu) Használati utasítás	(bg) Инструкции за употреба
(ca) MÓDE D'EMPLOI	(tr) KULLANIM TALimatları	(uk) УПАТСТВО ЗА УПОТРЕБА
(es) INSTRUCCIONES DE USO	(pl) INSTRUKCJA UŻYWANIA	(ja) 使用説明
(eus) INSTRUCCIÓN DE USO	(cs) NÁVOD K POUŽITÍ	(ko) 사용 지침
(pt) INSTRUÇÕES DE UTILIZAÇÃO	(sk) NÁVOD NA POUŽITIE	(zh) 使用說明
(dhr) INSTRUKÇÃO DE USO	(hr) UPUTE ZA UPORABU	(zh) 使用说明
(nl) GEBRUIKSAANWIJZING	(si) NAVODILA ZA UPORABO	(is) NOTKUNARLEIÐBEININGAR
(sv) BRUKSANVISNING	(sl) INŠTRUKCIJE DE UTILIZARE	
(fi) KÄYTTÖOHJEET	(et) KASUTUSJUHEND	

(en) ENGLISH

Indications: The products in this range are intended for use with a Colostomy, Ileostomy or Urostomy.
Intended Use: The product is intended for the management of stomal output.
Contraindications: None declared.
Precautions and Warnings:

- Check for sharp edges on the skin barrier before applying the pouch.
- Patients with a flush urostomy, when associated with megaureters and the absence of a urinary reservoir, should know that there have been reported cases of adhesive found in the urinary tract. Initiate appropriate clinical management and evaluate the continued use of the adhesive product if a mass is detected.
- This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection and cross-contamination. Physical properties of the device may no longer be optimal for intended use.
- A healthcare professional should assess the need for skin barriers with convexity prior to use.

Target Population: Ostomates requiring a convex baseplate.
Intended Users: All variants are intended for use by: Health Care Professionals, Caregivers and Ostomates.
Use Environment: All variants may be used in domestic and healthcare settings.
Incident reporting statement: If during the use of this device or as a result of its use a serious incident has occurred please report it to the manufacturer and to your national authority.
Disposal Statement: After use, this product may be a potential biohazard. Dispose of the pouch according to local disposal regulations.
© 2024 Convatec
™/® indicates a trademark of the Convatec Group

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