Hydrofiber™ Surgical Cover Dressing / Chirurgische Hydrofiber™ Wundauflage / Medicazione Hydrofiber™ Avanzata per la Ferita Chirurgica / Pansement Hydrofiber™ chirurgical de recouvrement / Adesivo de Cobertura Quirúrgico de Hydrofiber™ / Penso Cirúrgico de Hydrofiber™ / Hydrofiber™ Chirurgisch Dekverband / Hydrofiber™ förband för Postoperativ bandagering / Hydrofiber™ Postoperativtissidos / Hydrofiber™ postoperativ bandage / Hydrofiber™ bandage för postoperativ bandagering / Χειρουργικό επίθεσος Hydrofiber™ / غشاء Hydrofiber™

INSTRUCTIONS FOR USE / GEBRAUCHSANWEISUNG / ISTRUZIONI PER L'USO / NOTICE D'UTILISATION / INSTRUCCIONES DE USO / INSTRUCÇÕES DE USO / GEBRUIKSAANWIJZING / BRUKSANVISNING / KÄYTTÖOHJE / BRUGSANVISNING / BRUKSANVISNING / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / تعليمات الإستعمال

ENGLISH

PRODUCT DESCRIPTION
Aqualac™ Surgical cover dressing is a sterile post-operative dressing comprising of an inner (wound contact) non-woven pad composed of Hydrofiber™ technology stiller and absorbent. The pad is held in place behind the wound by a soft, top layer of Hydrocolloid which is bound to an outer polyurethane film and a windowed end contact layer consisting of polyurethane film sandwiched between two further layers of Hydrocolloid.

The outer polyurethane film of Aqualac™ Surgical cover dressing provides a bacterial and viral barrier to reduce the risk of infection. provided the dressing remains intact and there is no leakage.

This cover contains porcine gelatin.

INTENDED USE
Aqualac™ Surgical Cover Dressings have been designed to be used as a primary dressing.

Aqualac™ Surgical Cover Dressings may be used with the direction of a healthcare professional, for the management of acute, surgical wounds.

INTENDED USERS
Aqualac™ Surgical Cover Dressings are intended to be worn by patients who have undergone surgery.

TARGET PATIENT POPULATION
Aqualac™ Surgical Cover Dressings are intended to be worn by patients who have undergone surgery.

CLINICAL BENEFITS
Aqualac™ Surgical Cover Dressings provide a physical microbial, viral, and waterproof barrier to protect the wound.

Aqualac™ Surgical Cover Dressings absorb wound fluid and bacteria, providing a moist wound healing environment, whilst contouring to the wound/sutures/staples and removing dead-space between the wound and dressing interface.

Aqualac™ Surgical Cover Dressings are flexible and extensible to reduce the risk of dehiscence and aid skin movement.

INDICATIONS
Aqualac™ Surgical Cover Dressings are indicated for the management of acute, surgical wounds.

CONTRAINDICATIONS
Aqualac™ Surgical Cover Dressings should not be used on individuals who are sensitive to or who had an allergic reaction to the dressing or any of its components including colophony or its derivatives.

PRECAUTIONS AND WARNINGS
1. Sterility is guaranteed unless the pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or open to use and dispose of the device according to local regulations.

2. Aqualac™ Surgical Cover Dressings are for single use only and should not be re-used. Reuse may lead to increased risk of infection and cross-contamination and delayed healing.

3. This dressing should not be used with other wound care products without first consulting a healthcare professional.

4. Aqualac™ Surgical Cover Dressings are not compatible with oil-based products. Due to the sterilisation process there may be a slight odour on opening the primary packaging.

5. A healthcare professional should be consulted if any of the following are observed during dressing changes, irritation (reddening, inflammation), maceration (whitening of skin), hypergranulation (excesses of granulation), signs of infection (increased pain, bleeding, warmth/redness of surrounding tissue, pus, exudate), a change in wound colour and/or odour.

6. Aqualac™ Surgical Cover Dressings are designed to be MRH-Safe.

7. Aqualac™ Surgical Cover Dressings are not designed to be cut.

8. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local state and federal laws and regulations.

9. If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the competent authority of the Member State in which the user/s and/or patient is established.

10. Aqualac™ Surgical Cover Dressings are water resistant however, should not be submerged during bathing or swimming.

ASSOCIATED PRODUCT & THERAPIES
Aqualac™ Surgical Cover Dressings are not designed to be used in conjunction with any other products or therapies however, clinical judgement should be used to determine appropriate use where casts or splints are applied.

CONTACT AND DURATION OF USE
Aqualac™ Surgical Cover Dressings can be worn for up to 7 days, dressings should be changed earlier if clinically indicated. The Aqualac™ Surgical Cover Dressings are suitable for continuous wear until healing is achieved or an alternative wound management is clinically indicated.

DIRECTIONS FOR USE
To dress knee incisions, knee should be approximately at a 30-degree angle for optimal dressing application.

1. Before applying the dressing, cleanse the wound/incision site with an appropriate wound cleanser, rinse well with sterile saline and dry the surrounding skin.

2. Choose a dressing size to ensure direct contact between the wound/incision and the adhesive surface of the dressing.

3. Remove the dressing from the sterile pack.

4. Initially remove 3/4 of the backing paper from along the length of the dressing minimising finger contact with the central pad and the adhesive surface. Place at one end of the wound/incision line and gently roll the dressing over the wound/incision line and mould into place.

5. Remove the remaining 1/4 backing paper and gently mould the dressing over the wound/incision line in place for a secure adhesion to the skin.

6. The translucent hydrocolloid backing allows the clinician to monitor the central Hydrofiber™ technology pad and assess when the dressing needs to be changed. The dressing can be left in place for up to 7 days subject to regular clinical assessment and local dressing protocol.

7. All wounds/incisions should be monitored frequently. Remove Aqualac™ Surgical cover dressing when clinically indicated (i.e. leakage, excessive bleeding, suspicion of infection, etc).

Dressing Removal
To remove the dressing, press down on the skin with one hand and carefully lift an edge of the dressing with your other hand. Stretch the dressing to break the adhesive seal and remove.

Store at room temperature (10°C - 25°C/50°F - 77°F). Keep dry.

https://ec.europa.eu/tools/eudamed. This is the Summary of Safety and Clinical Performance (SSCP) location after the launch of European Database on Medical Devices/ EUDAMED.

If further information or guidance is needed, please contact Convatec Professional Services.

Made in Dominican Republic
© 2024 Convatec.
™/® indicates a trademark of the Convatec Group.

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