AQUACELTM Foam

INSTRUCTIONS FOR USE / ISTRUZIONI PER L'USO / NOTICE D'UTILISATION / INSTRUCCIONES DE USO / INSTRUCÇÕES PARA UTILIZAÇÃO / ГРУБКАСИВНИК / BRUKSANVISING / ОДНІЄНІЯ КАРТТООНЕЕТ / BRUGSANVISING / BRUKSANVISING / ОДНІЄНІЯ ХРХЗН / تعليمات الاستخدام

Adhesive and Non Adhesive AQUACELTM Foam Dressing / Medicazione HydrofoberTM adesiva e non adesiva in schiuma / Pansement hydrocolleux stérile issu de la technologie HydrofoberTM adhésif et non-adhésif / Απίσος de Espuma de HydrofoberTM Adesivo y No Adesivo / Pensó de espuma de Tecnologia HydrofoberTM aderente e não aderente / Adhesief en niet-adhesief de technologie HydrofoberTM schuimverband / Váhihtatande adheesive en niet-adheesive schuimverband met HydrofoberTM / Клейкий и неклеющий скаттооннеет с HydrofoberTM / Клейкий и неклеющий скаттооннеет с HydrofoberTM / Kiebande og ikke-kiebande skumbandage med HydrofoberTM / Teknologia / Kiebande og ikke-kiebande HydrofoberTM skumbandasje / Αυτοκόλλητο και μη-αυτοκόλλητο Απίσος με HydrofoberTM / رتبه‌تو لاصقه وغير لاصقه

ENGLISH

PRODUCT DESCRIPTION
AQUACELTM Foam dressing is a sterile, HydrofoberTM foam wound dressing, consisting of:
• a waterproof, breathable outer polyurethane film which provides a bacterial and viral barrier
• a soft, multi-layered, absorbent central pad containing polyurethane foam and HydrofoberTM (Sodium carboxymethylcellulose) wound contact layer. The pad absorbs and retains moisture, exudates and bacteria.
• AQUACELTM Foam Adhesive has a silicone adhesive skin contact border to adhere the dressing to the patient.

CLINICAL BENEFITS
AQUACELTM Foam is designed to protect fragile skin from moisture, shear and friction damage. AQUACELTM Foam is designed to provide a moist wound healing environment, aiding autolytic debridement and removing dead-space between the wound and dressing interface. AQUACELTM Foam is designed to manage excess exudate levels which may further damage the wound bed and surrounding skin.

INTENDED USE
AQUACELTM Foam dressings may be used, with the consultation of a healthcare professional, for the management of exuding, non-exuding wounds and protection of intact skin. AQUACELTM Foam dressing may be included in a comprehensive protocol of care to protect intact skin against breakdown.
AQUACELTM Foam dressings may be used as a primary or secondary dressing.

INDICATIONS
AQUACELTM Foam dressing is indicated for the treatment of:
• Leg ulcers
• Pressure Ulcers
• AQUACELTM Foam is not intended for use on patients with open, bleeding, or infected wounds.
• Partial thickness burns
• Traumatic wounds

INTENDED USERS
AQUACELTM Foam dressings are to be used by health care professionals, patients and carers.

TARGET PATIENT POPULATION
AQUACELTM Foam dressings are intended to be used on patients with one of the wound types listed in the indications or those who require skin protection from breakdown.

CONTRAINDICATIONS
AQUACELTM Foam dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or components stated above.
AQUACELTM Foam dressings are not compatible with oxidising agents such as hydrogen peroxide or hypochlorite solutions.

PRECAUTIONS AND OBSERVATIONS
Sterility is guaranteed unless the pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or opened prior to use and dispose of the device according to local regulations.
Do not use after stated expiry date.
AQUACELTM Foam dressings are for single use only, do not reuse. Reuse may lead to increased risk of infection and cross-contamination and delayed healing.
For highly exuding wounds, AQUACELTM Foam may be used as a secondary dressing in conjunction with a primary dressing, such as AQUACELTM Extra.
This dressing should not be used with other wound care products without first consulting a healthcare professional.
If the adhesive dressing is cut, additional tape may be required to secure it in place.
During the use of this device or as a result of use a serious incident has occurred, please report it to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

CONTACT AND DURATION OF USE
AQUACELTM Foam dressings can be worn for up to 7 days, dressings should be changed earlier if clinically indicated, and are safe for continuous use.

DIRECTIONS FOR USE
1. **Storage**
Store at room temperature (10°C - 25°C/50°F - 77°F). Protect from light. Keep dry.
2. **Wound Site Preparation and Skin Cleaning:**
Before applying the dressing, cleanse the area with an appropriate wound cleanser and dry the surrounding skin.
3. **Dressing Preparation and Application:**
a. Before applying to the wound, ensure the dressing size to ensure that the central absorbent pad is at least 1 cm (1/3 inch) larger than the wound.
b. When applying to intact skin choose a dressing which covers the entire vulnerable area.
c. Remove the dressing from the sterile pack.
d. AQUACELTM Foam Adhesive only: Remove the release liner, minimising finger contact with the adhesive surface.
e. Hold the dressing over the wound or area of skin to be protected and lift the centre of the dressing with the centre of the wound or area of skin to be protected.
f. Ensure the dressing is not stretched upon application. For AQUACELTM Foam Adhesive smooth down the adhesive border. Ensure joints are sufficiently flexed to allow patient mobility.
g. The dressing may be cut; if cut, a residual tape or film dressing should be used to create a seal and ensure a bacterial barrier.

4. **Dressing Removal:**
a. The dressing should be changed when clinically indicated (the dressing should not be allowed to become saturated or if it becomes damaged/soiled).
b. If removing the dressing is difficult, the dressing should be saturated with sterile saline or water and removed slowly.
c. To remove the dressing, peel down gently on the skin and carefully lift one corner of the dressing, continue until all edges are free. Carefully lift away the dressing and discard as per local clinical protocols.
d. Only during preventative care, dressing corner can be lifted for inspection of skin and reassessed. Dressing must be changed if skin is broken.

Visit us online at www.convatec.com.
Made in UK
© 2023 Convatec Inc.
**Indicates a trademark of Convatec Inc. AQUACEL, the AQUACEL logo and Hydrofober are registered trademarks

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