Kaltostat™

Calcium-Sodium Alginate Wound Dressing / Calcium-Natriumalginat Wundverband / Medicazione Per Ferite in Alginato Di Calcio E Sodio / Pansement d'Alginate Calcium-Sodium / Apósito de Alginato Cálcio-Sódico Para el Tratamiento de Heridas / Apósito de Alginato de Cálcio-Sódio para Tratamento de Feridas / Calcium-Natrium Alginat Wundverband / Calciumnatriumalginat Sårförband / Kalsiumnatriumalginattihaavavädos / Calcium-Sodium Alginatne Sårbandage / Calcium-natrium Alginat Sårbandasje / Alyuved Enfüqes Tenqumayur Apqektow / Netprow / مضاد الكالسيوم صوديوم الجيلاتن المرحوح الماصحة

INSTRUCTIONS FOR USE / GEBRAUCHSANWEISUNG / ISTRUZIONI PER L'USO / NOTICE D'UTILISATION / INSTRUCCIONES DE USO / INSTRUÇÕES DE UTILIZAÇÃO / GEBRUIKSAANWIZING / BRUKSANVISNING / KÄYTTÖOHJEET / BRUGSANVISNING / BRUKSANVISNING / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / تعليمات الاستخدام

PRODUCT DESCRIPTION
Kaltostat™ Wound and Cavity Dressings are soft, white-to-off-white, sterile dressings of calcium-sodium alginate fibre. Kaltostat™ is available as a non-woven flat wound dressing and as a 2 gram wound cavity dressing shaped like a 'rope' (Kaltostat™ Rope).

INTENDED USE
Kaltostat™ dressings may be used, with consultation of a healthcare professional, for the management of moderate to heavily exuding wounds. Kaltostat™ dressings are designed to absorb exudate and protect the wound from contamination.

Kaltostat™ Rope dressing is designed to be used to manage cavity wounds.

INTENDED USERS
Kaltostat™ dressings are intended to be used by healthcare professionals who are suitably qualified and trained in their use.

TARGET PATIENT POPULATION
Kaltostat™ dressings are indicated for use in adults and children over 12 months of age.

CLINICAL BENEFITS
Kaltostat™ dressings absorb exudate and promote haemostasis. Kaltostat™ dressings contain calcium and the availability of calcium ions at an injury site helps to promote haemostasis.

Kaltostat™ dressings absorb fluid directly into the fibres. Exudate is absorbed into the dressing on contact with the wound, keeping fluid away from the peripheral skin. The risk of wound edge and surrounding tissue maceration is therefore reduced.

Kaltostat™ dressings are also designed for removal with little or no damage to the newly formed tissue.

Kaltostat™ Rope can be used to manage cavity and deep wounds.

INDICATIONS
Kaltostat™ dressings are indicated for the management of the moderate to heavily exuding wounds listed below:

- Pressure ulcers, venous leg ulcers, arterial leg ulcers, mixed aetiology ulcers
- Diabetic foot ulcers
- Surgical wounds, traumatic wounds, bleeding wounds
- Skin graft donor sites

CONTRAINDICATIONS
Kaltostat™ dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components including calcium sodium alginate or other seaweed derivatives.

Kaltostat™ dressings are not indicated for dry wounds or third-degree burns. Kaltostat™ dressings are not intended to control heavy bleeding. Alternative measures must be considered in those emergency situations where large quantities of blood may be lost.

PRECAUTIONS AND OBSERVATIONS

- The use of Kaltostat™ dressings is not recommended for infants less than 12 months of age due to the presence of manganese within the dressing (s 20 ppm).
- Kaltostat™ dressings are not indicated for use as a surgical sponge.
- Kaltostat™ dressings are single use devices and should not be reused. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use.
- Stability is guaranteed unless the pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or open prior to use and dispose of the device according to local regulations.
- Kaltostat™ dressings should not be used with other wound care products without first consulting a healthcare professional.
- Should you observe irritation (reddening, inflammation), maceration (softening of skin), hypergranulation (excess tissue formation) or sensitivity (allergic reaction), consult a healthcare professional.
- Where Kaltostat™ dressings have initially formed a gel that is allowed to dry out, rehydration of the dressing with sterile saline will be required to enable painless removal.

The dressing may be used on infected wounds only under the care of a healthcare professional.

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local state and federal laws and regulations.

If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

ASSOCIATED PRODUCT & THERAPIES
Kaltostat™ dressings require fixation with an appropriate secondary dressing and can be used in conjunction with compression therapy.

CONTACT AND DURATION OF USE
Kaltostat™ dressings can be worn for up to 7 days, dressings should be changed earlier if clinically indicated.

DIRECTIONS FOR USE

1. Wound Site Preparation and Cleansing
Before applying Kaltostat™ dressings, cleanse the areas with an appropriate wound cleanser.

Under the care of a healthcare professional, if necessary, before applying Kaltostat™ dressings, the wound should be debrided of excessive necrotic tissue and eschar and irrigated with an appropriate wound cleanser.

2a. Dressing Preparation and Application (Kaltostat™ Wound Dressing)

- Kaltostat™ dressings should be trimmed to the exact size of the wound avoiding any overlap onto the surrounding skin.
- Kaltostat™ dressings should be applied dry onto or into the wound. An appropriate secondary dressing should be used to secure Kaltostat™ dressings in place.
- Discard any unused portion of the product after dressing the wound.

2b. Dressing Preparation and Application (Kaltostat™ Rope Only)

- Kaltostat™ Rope dressing can be cut to the required length to fill the cavity wound. Do not pack too tightly.
- Kaltostat™ Rope dressing should be applied dry onto or into the wound. An appropriate secondary dressing should be used to secure Kaltostat™ dressings in place.
- Discard any unused portion of the product after dressing the wound.

3. Dressing Change and Removal

- Change the Kaltostat™ dressings when strike-through of the secondary dressing occurs or as indicated by good clinical practice.
- Removal should be trouble free on heavily exuding wounds. Kaltostat™ dressings will gel at the wound/dressing interface and lift away.
- As with any dressing, Kaltostat™ dressings should be removed from the wound and the wound cleansed at appropriate intervals in accordance with good clinical practice, or up to a maximum of 7 days.
- After removal, dressings should be discarded as per local clinical protocols.

Store in cool dry place. (10°C/25°C/50°-77°F). Keep Dry. Protect from light. If further information or guidance is needed, please contact Convatec Professional Services.

™/® Indicates a trademark of the Convatec Group

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© 2025 Convatec