INSTRUCTIONS FOR USE / GEBRUIKSANWEISUNG / ISTRUZIONI PER L'USO / NOTICE D'UTILISATION / INSTRUCCIONES DE USO / INSTRUÇÕES PARA UTILIZAÇÃO / GEbruIKSINSTRUCTIES / BRUKSANVISNING / KÄYTTÖOHJEET / BRUGSANVISNING / BRUKSANVISNING / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

Translations available: [Translations](#)

Adhesive and Non Adhesive Hydrofiber™ Foam Dressing with Silver / **AQUACEL™** Ag Schaumverband mit Silber / **Medicazione Hydrofiber™** adesiva e non adesiva in schiuma con ion di argento / **Pansement hydrocellulaire à l'argent**, issu de la technologie **Hydrofiber™** adhésif et non-adhésif / **Apósito de espuma de Hydrofiber™** Adhesivo y no Adhesivo / **Ag Schaumverband mit Silber / Medicazione Hydrofiber™** adesiva e non adesiva em schiuma com iónio de prata / **Hydrofiber™** teknologi og salv / **Klebefolie und ikke-klebefolie Hydrofiber™** skumbandsje med silver / **Autoklātētā kārtīgā un nēkārtīgā Hydrofiber™** ēriņš / **Επίθεση Αργύρου με Hydrofiber™** / **مطاطة ملتحمة وملتحمة غير ملتحمة حبر فلز** / **Hydrofiber™** باعجية قصبة وقصبة غير قصبة (Hydrofiber™)

(EN) ENGLISH

PRODUCT DESCRIPTION

AQUACEL™ Ag foam dressing, adhesive and non adhesive, is a range of sterile Hydrofiber™ foam wound dressings consisting of a waterproof outer polyurethane film and a multi-layered absorbent pad with the adhesive dressing having a silicone adhesive border. The multi-layered absorbent pad contains a layer of polyurethane foam and a non-woven wound contact layer of Hydrofiber™ technology (sodium carboxymethylcellulose). The dressing's Hydrofiber™ wound contact layer contains 1.2% w/w ionic silver.

The silver in the dressing kills a broad spectrum of wound bacteria, and helps reduce the risk of wound infection. The outer film layer provides a waterproof viral and bacterial barrier which protects the wound from external contaminants, reducing the risk of infection. The film also helps to manage the moisture vapor transmission of the exudate absorbed by the dressing.

The foam and Hydrofiber™ materials within the pad absorb high amounts of wound fluid and bacteria. The Hydrofiber™ wound contact layer creates a soft cohesive gel that intimately conforms to the wound surface, maintains a moist environment in the wound (which supports the body's healing process) and aids in the removal of non-viable tissue from the wound (autolytic debridement) without damaging tissue. The dressing's non-adherent wound contact layer allows for minimal trauma on removal. The adhesive dressings in the range have a silicone border which provides secure, skin-friendly adhesion and supports non-traumatic dressing removal.

AQUACEL™ Ag foam dressing may be used as a primary or as a secondary dressing. It may be used alone or in combination with other wound care products and can be cut to shape and size, as directed by a healthcare professional. If the adhesive dressing is cut, additional tape may be required to secure it in place. Additional tape, or other means of fixation, will be required to secure the AQUACEL™ Ag foam non adhesive dressing in place.

AQUACEL™ Ag foam adhesive dressing's film backing acts as a barrier to the wound against bacterial and blood-borne viral pathogens (e.g. HIV and Hepatitis viruses). The use of this dressing neither guarantees nor warranties against AIDS or Hepatitis virus transmission.

INDICATIONS

Under medical supervision AQUACEL™ Ag foam dressing range may be used as stated in the DIRECTIONS FOR USE for the management of both chronic and acute wounds such as:

- wounds where there is an infection or an increased risk of infection;
- partial thickness (second degree) burns;
- diabetic foot ulcers, leg ulcers (venous stasis ulcers, arterial ulcers and leg ulcers of mixed etiology) and pressure ulcers/sores (partial & full thickness);
- surgical wounds;
- traumatic wounds;
- wounds that are prone to bleeding, such as wounds that have been mechanically or surgically debrided;
- oncology wounds with exudate, such as fungoelosive tumors, fungating carcinoma, cutaneous metastasis, Kaposi's sarcoma, and angiiosarcoma.

CONTRAINDICATIONS

AQUACEL™ Ag foam dressings should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.

PRECAUTIONS AND OBSERVATIONS

Caution: Sterility is guaranteed unless pouch is damaged or opened prior to use.

This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use.

Clinicians/Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.

During the body's normal healing process, non-viable tissue may be removed (autolytic debridement) from the wound, which could make the wound appear larger after the first few dressing changes.

Should you observe irritation (reddening, inflammation), maceration (whitening of skin), or hypergranulation (excess tissue formation) consult a healthcare professional.

The wound should be inspected during dressing changes for (1) signs of infection (increased pain, bleeding, warmth/fredness of surrounding tissue, wound exudate), (2) there is a change in wound colour and/or odour, (3) any other unexpected symptoms occur (e.g. maceration, or hypergranulation).

Appropriate supportive measures should be taken where indicated (e.g. use of graduated compression bandaging in the management of venous leg ulcers or pressure relief measures in the management of pressure ulcers). The control of blood glucose, as well as appropriate supportive measures, should be provided to persons with diabetic foot ulcers.

AQUACEL™ Ag foam dressing, adhesive and non adhesive, was determined to be MR-Safe according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

DIRECTIONS FOR USE

1. If the immediate sterile product pouch is damaged, do not use.

2. Wound Site Preparation and Cleansing:

Before applying the dressing cleanse the wound area with an appropriate wound cleanser and dry the surrounding skin.

3. Dressing Preparation and Application:

a. Choose a dressing size and shape to ensure that the central absorbent pad (area within the adhesive window), is 1cm (1/3 inch) larger than the wound area.

b. Remove the dressing from the sterile pack, minimizing finger contact with the wound contact surface and the adhesive surface where applicable. Remove the release liner if using the adhesive dressing.

c. The dressing can be cut to shape for convenience.

d. Hold the dressing over the wound and line up the center of the dressing with the center of the wound. Place the pad directly over the wound. For the adhesive dressing, smooth down the adhesive border.

e. An appropriate retention bandage or tape should be used to secure the dressing in place if the dressing does not have an adhesive border, or if the adhesive dressing has been cut.

f. For difficult to dress anatomical locations, such as the heel or the sacrum, the specially shaped adhesive dressings may be used.

g. Discard any unused portion of the product after dressing the wound.

4. Dressing Removal:

a. The dressing should be changed when clinically indicated (e.g. leakage, bleeding, increased pain). Maximum recommended wear time is up to seven days.

b. The wound should be cleansed at appropriate intervals.

c. To remove the dressing, press gently on the skin and carefully lift one corner of the dressing. Continue until all edges are free. Carefully lift away the dressing and discard according to local clinical protocols.

Store at room temperature (10° C - 25° C/50° F - 77° F). Keep dry. Protect from light.

If further information or guidance is needed, please contact ConvaTec Professional Services or visit us online at www.convatec.com.

Made in UK

© 2021 ConvaTec Inc.

*™ indicates a trademark of ConvaTec Inc.

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