

ARTWORK SPECIFICATION

Artcode	1730227V1
Artwork Description	INS AQUACEL AG+ EXTRA (NON-DE) EUR
ICC*	413567
SAP*	1708331
KOM/REF*	N/A
IS Code*	N/A
Component Type	Insert
LR#	LR-004775
Existing Art Code being Revised	1730227A1

REQUIREMENT	SPECIFICATION
Bucket	EUR
Local Code(s)	N/A
Size Code / Dieline	N/A
Dimensions	445mm x 810mm FLAT, Folds: Height – 90mm Width – 90mm
Color(s) / Colour(s)	PMS 527
Specification Reference(s)	MT56-001
Bar Code: UPC (A) or GTIN**	N/A
Bar Code: PZN – Code 39**	N/A
Bar Code: JAN-13**	N/A

* NOTE: Bar Code: See artwork for positioning, magnification and truncation. Scanability verified to ANSI grade "C" or better.

HISTORY

14-Sep-2020 / 1730227 : New insert created based on existing art 1708881G1 with the following changes: German language removed, color changed from Black to blue (PMS647), CE numb.updated from CE0086 to CE2797. All symbols updated per ISO15223-1 and TM added after Hydrofiber tradenames. CVT Logo and AQUACEL updated. Authorized representative in the European Community symbol with address added after DEE address.Representant of Legal (Clare Bates) was informed about change via email 14-Sep-2020.

20-Dec-2021 / 1730227A1: 1730227A1: Burn indication added. Symbol "Importer" added next to EC REP address. Font size can be shrink no less than 6pt. Symbols Manufacturer, Caution, Keep away from sunlight, Does not contain natural rubber latex added and explained. Explanation of these symbols "REF", "LOT", "Use-by date", "Do not re-use", "Consult instruction for use", "Keep dry", "Temperature limit", "Sterilized using irradiation", and "Do not use if package is damaged" updated acc. to SM-079379. Versiva™ XC™ and shortcut (BeCl) removed.Representant of Legal (Gregg Vicinanza) was informed about these changes by email (12-Jan-2022).

21-Nov-2022 / 1730227V1: All current symbols and their explanations updated. Font size shrank to 6pt. Version updated to V1. Caution symbol with explanations removed. Language titles updated. ConvaTec and AQUACEL Ag Extra logo updated according Master Brand req. Addresses and phone numbers updated. Texts were revised. CH REP symbol and address with Importer symbol added. "UKCA 0086" added. Symbols with their explanations in all languages added: "MD", "MR - Safe", "EC REP", "Date of manufacture", "UDI", "Contains a medicinal substance", "Single sterile barrier system", "Single sterile barrier system with protective packaging outside", "CH REP", "Importer". Color changed to PMS 527.